



# PACS Pro™ Standard RSVS Viewer

Unlike most viewers on today's market, RSVS Viewer is Diagnostic-quality viewing software that can support up to 5MP Medical Grade LCDs. RSVS can be used as stand-alone viewing software on any workstation that requires a DICOM viewer. A powerful, customizable, easy-to-use viewer that allows users to configure the toolbar and preset mouse buttons to their preferences, RSVS allows you to easily access the tools you use most often. For use on all review stations and Radiologist workstations.



## RSVS Viewer Highlights:

- Scout/Reference Lines
- Multi-Planner-Reconstruction
- KIN (Key Image Note)
- GSPS/Bitmap Overlays
- Automatic Alignment
- Study Comparison (New/Old)
- Playback (Multi-Frame)
- DICOM Scan/Send

## Patient Worklist Features:

- Patient/Study/Series Levels
- Image Sorting
- Search Patient Worklist
- Single/Multiple list Selection
- Adjacent/Non-Adjacent Images
- Individual Image Selection
- Image Storage
- House Keeping Utilities

## Measurement Tools:

- Distance
- Angle
- Area
- Density
- Histogram
- Erase Function
- Calibrate
- Measure Heart/Chest Ratio

## DICOM Network:

- Peer to Peer Transfer
- Image Transfer
- Query Retrieve
- DICOM Network Manager
- DICOM Communication Activities

## Viewing Features:

- Zoom In/Out
- Window Level
- Pan/Roam
- Pseudo Color
- Image Reset
- Display Multiple Studies
- Easy Configuration Utility
- JPEG Lossless/Lossy
- Auto Image Routing
- Color Display Support
- Event Log
- Mapping Algorithm Selection
- Magnify/Cine/Shutter/Crop
- Export Image
- Print (Paper/Film)
- Information Text Overlay
- Filters
- Preset
- Import/Display image from CD
- Diagnostic Embedded Reports
- Study Notes
- Annotations
- Attach document to study
- Image orientation
- Cross link
- Mammography display layout

