



IDD X-Ray Film and Paper Digitizer S6000 Series

Innovative Digital Devices (IDD) S6000 film and paper digitizer can be used in a number of user environments. The applications include but are not limited to PACS, Teleradiology, Orthopedic surgery, Podiatry and Veterinary practices. S6000 digitizers are the only medical digitizers capable of processing both film and paper, eliminating the need for an additional paper document scanner.

IDD digitizers offer more “up time” and afford easy maintenance to enhance user productivity.

IDD S6000 Series Specifications

600 dots per inch resolution (up to 8K by 10K)

Scan many different film and paper sizes, up to 14x35”

8-bit/16-bit processing (65536 shades of gray)

Gamma, brightness, contrast, and exposure controls

Optical Density of better than 4.0

Light and compact, only 19.8 lbs.

Wall mountable

Feed tray and catch basket

Optional Multiple Film Feeder

Pictured without Film Auto-Feeder



Pictured with Optional 15-film Feeder



Increase Productivity with Multiple Film Feeder:

- Digitize multiple X-Ray films with one click
- 10 to 15 mixed-size films
- 50-60 films per hour