

KODAK HYPER SPEED G MEDICAL FILM



Cut patient exposure up to 50% without sacrificing image quality or changing your workflow

KODAK Hyper Speed G Medical Film allows patient and staff radiation exposure to be reduced by up to one-half while maintaining high diagnostic image quality. This unique high-speed general-purpose film features zero-crossover technology for high resolution and increased visibility of fine detail. Its speed reduces image blurring and the need for retakes due to patient motion.



IMAGING APPLICATIONS

This film's low-dose advantage is appropriate for general purpose imaging including:

- Pediatric examinations*
- Obese patient examinations
- Spine/scoliosis examinations
- Portable examinations
- Chiropractic examinations

**For infants, KODAK INSIGHT Pediatric Imaging Film is recommended.*

HIGH SPEED AND HIGH QUALITY

Hyper Speed G Medical Film delivers twice the speed of premier KODAK X-SIGHT G/RA Film—with equivalent sharpness and similar high-contrast curve.

KODAK MEDICAL SYSTEMS



KODAK HYPER SPEED G MEDICAL FILM

USE WITH YOUR GREEN SCREENS

Hyper Speed G Medical Film can be used with any green-light-emitting screens. There is no need to invest in screens or cassettes.



WE HEAR YOU!


KODAK Hyper Speed G Medical Film is the first high-resolution radiographic film in 20 years to offer significant reduction in patient dose. Developed in direct response to customer needs and requests, it is another milestone in Kodak's continued investment in radiographic film systems that set new standards for performance and image quality.

MAINTAIN YOUR WORKFLOW


Conversion to this high-speed film requires no change in workflow—so you maintain efficiency and productivity.

MORE INFORMATION

Learn more about KODAK Hyper Speed G Medical Film—and other KODAK radiographic film offerings. Contact your local Kodak representative or visit us on the Web at www.kodak.com/go/medical.



- Hyper Speed G/X-SIGHT Screen
- Hyper Speed G/Regular Screen
- Hyper Speed G/Medium Screen
- T-MAT H/RA/Regular Screen
- T-MAT G/RA/Regular Screen
- X-SIGHT G/RA/X-SIGHT Screen



Health Group
EASTMAN KODAK COMPANY
343 State Street
Rochester, NY 14650
1-877-TO-KODAK (1-877-865-6325)

KODAK CANADA INC.
3500 Eglinton Avenue West
Toronto, Ontario M6M 1V3
CANADA

Outside the U.S. or Canada, please contact your local Kodak company.

www.kodak.com/go/medical

HEALTH GROUP
A BETTER VIEW OF LIFE.

