

MIN-R Mammography Processor

The MIN-R mammography processor delivers quality processing in an affordable, compact package. A specially designed film path and dryer improve image quality by substantially reducing the incidence of processing and drying artifacts. The result is increased diagnostic confidence with fewer repeats and more time with patients.



Other Features:

- Automatic standby control turns off the film drive, dryer blower, and dryer heater to save energy when no film is being processed
- Digital developer temperature readout simplifies regulatory record keeping
- The 142-second processing cycle (standard cycle) delivers 164 mammography film sheets per hour (18 x 24 cm)
- Drop-in" successor to a M35M processor; shares the same small footprint and electrical requirements, processor stand and replenisher tanks

Available Options:

- Stand
- Exhaust Fan Kit