



## PACS Pro™ Image Acquisition & Viewing Software


Intuitive software to easily enter Patient demographics, capture, view & distribute DICOM images from Carestream PoC, Vita & Orex ACL CR Scanners

**PACS Pro™ Image Acquisition & Viewing software is the optimal filmless solution for practices needing to: Acquire, View, Distribute and Manage DICOM images**


PACS Pro™ software is available exclusively through TI-BA Enterprises and its dealer partners. This software integrates with Kodak Point-of-Care (PoC), Vita and Orex ACL branded CR systems. PACS Pro provides an intuitive GUI making customer applications training simple and fast. Our “One-Click” design provides rapid access to many key functions with a single mouse-click; i.e. CD Burning, DICOM send, Patient Edit, DICOM Viewer, etc...

Acquisition software provides short-term image archive and viewing capability as well as the ability to send images over a DICOM network or distribute/archive via CD & DVD. PACS Pro™ Image Acquisition includes a built-in diagnostic quality DICOM image viewer with standard measuring and QC tools. It also supports an optional DICOM Modality Worklist (MWL) to register and retrieve patient demographics from your RIS/HIS, which saves time and reduces the chance of errors when keying patient information.


**Mobile Applications:** PACS Pro provides the ability to add up-to-eight unique DICOM tags as well the ability to attach color pictures and electronic documents to the study for future reference.



PACS Pro Patient Worklist Screen



PACS Pro Configuration Screen



PACS Pro DICOM Viewer with Measurement Toolbar Displayed

For Sales Contact Your X-ray Service Provider: