



## Rebuilt CR Systems

A complete CR image acquisition solution for lower volume medical practices seeking to make the digital transition.

Since 1979, TI-BA Enterprises has been the benchmark of quality, consistency and customer support in the refurbished medical imaging equipment industry. This celebrated tradition of workmanship and customer service has been assigned to the ever-evolving digital market.

Each TI-BA refurbished CR is fully inspected and run through a battery of tests before it is resold. These systems are equipped with new Dell PC's, loaded with TI-BA's PACS Pro Acquire software, tested and pre-staged for your convenience. All cassettes and screens are cleaned and inspected to ensure they meet our quality standards.

The end product is a true 'plug and play' CR System.



**Kodak CR Description:** The Orex ACL4 and its successor, the Kodak Point-of-Care 140, are ideal for private specialty practice, Veterinary, Imaging Centers, mobile installations and satellite locations. Durable, compact, and lightweight, these 40 Plate Per Hour (pph) CR readers enable all imaging functions to be performed efficiently at the point-of-patient-care with one affordable package price.

System features a single-plate CR reader, One (1) 10x12, One (1) 14x17 cassettes with flexible phosphor screens, advanced PACS Pro Image Acquisition & DICOM Viewing software and a PC-based review station with 19" LCD.

### Refurbished ACL4/PoC 140 CR System Specifications

<b>Throughput (cassettes per hour)</b>	Up to 41 - 14x17" Cassettes
<b>Grayscale Resolution Acquisition Display</b>	12 bits per pixel 12 bits per pixel
<b>Software</b>	Integrated PACS Pro™ Acquisition Software with up to 500GB Storage capability. DICOM: View, Send Store, QC, AutoSend, Patient Worklist, Multi-Patient CD/DVD Burn, DICOM Tags, Touch Screen Capable, DICOM 3.0 standard  TI-BA can offer Human & Veterinary Systems
<b>Dimensions (W x D x H)</b>	29" X 26" X 14"
<b>Weight</b>	88 lbs
<b>Dell PC Workstation Specs</b>	OptiPlex MiniTower, 2 GB RAM, 320GB HDD 19" LCD, DVD+/-RW, Windows XP Pro
<b>Power Requirements</b>	Single Phase 50-60 HZ, 200 VA, 100 AVC – 240 AVC (±10%), UPS Required
<b>Cassettes/IP's</b>	One (1) 10x12 Cassette w/Phosphor Plates One (1) 14x17 Cassette w/Phosphor Plates
<b>Warranty</b>	6 months CR & Cassettes – Parts 1 Year Software Support & Upgrades Dell 3 Year Warranty on PC



*Pictured with Optional TI-BA CR Cart*