



TI-BA Enterprises, Inc. 2011 Catalog

*TI-BA Enterprises, Inc., continues to offer the best in both analog and digital **refurbished and new** equipment, OEM-quality replacement and original parts and equipment training courses.*



Our full-line of *refurbished* Carestream Film Processors, Multiloaders, Orex ACL2's, ACL4's, Digitizers and OEM-quality replacement parts has positioned us as the leader in the marketplace for 30 years. We are committed to supplying both you and your customer with the very best in *refurbished and new* equipment. Please call for specifics and pricing!

As we approach our 31st year in the business, we remain committed to quality and complete customer satisfaction. Our technical support team is always ready to assist you, whether it is to identify a part, troubleshoot a problem, or hunt down and find an obsolete item.



As technology continues its rapid change, we've expanded our product line, offering our customers a wide and diverse range of systems.



From film processors, to digitizers to CR Systems, we offer a multitude of product lines and manufacturers, *affordably priced*, to help you and your customer transition into the new technology.





TI-BA Enterprises, Inc.

25 Hytec Circle
Rochester, New York 14606
585-247-1212
www.ti-ba.com



FI Sales, LLC



Table of Contents

LIST PRICES EFFECTIVE MARCH 1, 2011.....	
IMPORTANT STUFF!	IV
EQUIPMENT WARRANTY	V
PARTS WARRANTY	VI & VII
FILM PROCESSOR PARTS	8-41
MULTILOADER ACCESSORIES	44-43
DRY LASER ACCESSORIES	47-45
RACK KITS	46
SUPER RACK KITS	47
COMPLETE TRANSPORT SYSTEMS	48
TUBE KITS	49
TANK FITTING KITS	46
TRANSPORT SYSTEMS	51-52
PUMP AND MOTOR CHART	53
ORDER FORM.....	56
POC PARTS	58-61
GUIDE TO GUIDESHOES.....	62

IMPORTANT STUFF!

- 1) All “T” prefixed parts are independent of original manufacturer.
- 2) Part numbers prefixed with an “R” **may** contain *reconditioned part(s)*.
- 3) Part numbers with a **red asterisk (*)** are followed by their new, orderable substitute part number (**in red!**).
- 4) List prices shown are those suggested by the original manufacturer and are subject to change without notice.
- 5) Shipping is FOB, Rochester, New York (*unless otherwise specified*).
- 6) **Hours of Operation:** Monday thru Friday 7:30 a.m. – 5:30 p.m. **(EST)**
- 7) We have on-site Training Classes. Visit us online at: **<http://www.ti-ba.com>** for details!
- 8) There are no returns on publications. Once you buy it, it’s yours!
- 9) We are an **authorized Agfa, Carestream, Konica-Minolta and AFP dealer**.
- 10) We are the U.S. distributor for MAVIG, Jamieson and Vanguard CINE parts.

REFURBISHED EQUIPMENT WARRANTY

Terms and Conditions

TI-BA Enterprises' *refurbished equipment* is covered by a one-year *limited warranty* on parts only from date of invoice. TI-BA Enterprises will, *at its discretion*, replace or repair the product at no charge. Choice of a Unit Service Exchange will be solely at the discretion of TI-BA Enterprises, and will be based on communications between TI-BA Enterprises, Inc., and the end-user.

This warranty assures the equipment to be *free from manufacturing defects*★. *Any qualified warranty part* will be replaced at no charge.

★★★ Any equipment damaged in shipment or transit is not covered under warranty. If this situation occurs, YOU ARE RESPONSIBLE for submitting a claim to the carrier. Our equipment is labeled such that the carton should be opened and equipment inspected PRIOR TO SIGNING OFF on the shipment (while the delivery driver is still there). It is advisable before accepting any package/shipment that you open and inspect for concealed damage! ★★★

This **warranty excludes:**

- A. Normal wear and tear
- B. Lack of proper maintenance
- C. Improper installation or change in voltage after installation
- D. Cosmetic damage after installation
- E. Improper de-installation of part (*cut wires on motors, pumps, etc.*)

IN NO EVENT shall TI-BA Enterprises be liable for consequential or incidental damages. **This warranty excludes any and all "Acts of God" (or any of his "helpers").**

★To obtain warranty part(s), dealer must provide TI-BA Enterprises with:

- 1) The **TI-BA assigned serial number** of the unit and a new purchase order number for the replacement part(s). Replacement part(s) will be shipped as requested and will be billed at normal prices.
- 2) **Defective parts must be returned within 20 days** to TI-BA Enterprises in Rochester, NY, to qualify for replacement credit. Items not returned within 30-day period will be a chargeable item. No credit will be issued.
- 3) **Only** if part is found to be **defective in manufacture** upon inspection by TI-BA Enterprises, Inc., will credit be issued.

If part(s) are not returned within the 20-day limit, full payment of the invoice is expected.

PARTS WARRANTY

(Revised 01.08)

All parts orders are subject to acceptance by TI-BA's Customer Service Department in Rochester, NY.

PRICING

Prices are FOB Rochester, NY, unless specified otherwise, and are subject to change without notice. Prices do not include labor, installation or freight. There is no minimum order requirement, nor is there ever a "rush" or "expedite" charge.

For products or equipment contained within this catalog without pricing, please call for price quotation.

TERMS OF PAYMENT

Subject to Credit Department approval; 2% NET 10 Days; NET 20 Days from date of invoice.

PARTS WARRANTY

All items, with the exception of consumable items, supplied by TI-BA Enterprises, are warranted by TI-BA Enterprises to be **free of manufacturing defects and to function properly for three months from TI-BA's billing date**. TI-BA MAKES NO OTHER WARRANTIES, EXPRESS, IMPLIED OR OF MERCHANTABILITY FOR THESE ITEMS.

Repair and the option of replacement, in the event that TI-BA is unable to repair the item, are TI-BA's only obligations. TI-BA WILL NOT BE RESPONSIBLE FOR ANY CONSEQUENTIAL OR INCIDENTAL DAMAGES RESULTING FROM THE SALE OR USE OF THESE ITEMS, EVEN IF LOSS OR DAMAGE IS CAUSED BY NEGLIGENCE OR OTHER FAULT OF TI-BA.

All items sold by TI-BA Enterprises, Inc., are supplied without return privileges, except for items which are

returned for warranty service consideration, damaged in transit or shipped in error by TI-BA. Any other return requires authorization by our Customer Service Department and **if authorized**, are subject to inspection, handling and restocking charges. Requests for return **must include** part number being returned, description, original purchase order number (associated invoice) and reason for return.

Parts Warranty *(cont'd)*

PARTS ORDERED INCORRECTLY: May be returned within 20 days of the original order date. **Parts returned after 20 days will not receive credit.**

PARTS SHIPPED INCORRECTLY: May be returned within 10 days of the original order date. These parts are not subject to TI-BA restocking fees.

PARTS ORDERED FOR WARRANTY REPAIR: The defective parts must be returned within 20 days to receive warranty credit.

NO CREDIT WILL BE ISSUED for warranty parts returned after 20 days. THERE ARE NO EXCEPTIONS!!!!!!

ELECTRONIC PARTS: TI-BA Enterprises WILL NOT ACCEPT any return of electronic parts, if the original package is open and/or if the original package seal is broken.

NOTE: To facilitate the return of in-warranty TI-BA parts, please include a complete description of the reason(s) for the return. Items being returned for warranty consideration **MUST BE RETURNED within 20 days.**

Items returned after 20 days will void all return privileges to customer and any warranty obligations by TI-BA Enterprises.

All items should be returned, transportation prepaid by customer, to:

Customer Service Department

TI-BA Enterprises, Inc.
25 Hytec Circle
Rochester, New York 14606



Film Processor Parts

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
TL1481	TTL1481	Alignment Tool	\$11.74
TL2094	TTL2094	Red Test Lead	\$52.61
TL2095	TTL2095	Black Test Lead	\$52.61
TL2170	TTL2170	Pinch Clamp	\$253.45
TL2431	TTL2431	Air Meter w/case	\$190.81
TL4391	TTL4391	Cable-Laptop Interface	\$158.70
TL4398	TTL4398	Sm Chip Removal Tool	\$130.10
004435K	T004435K	Pump Housing Kit	\$153.09
004488TL	TTL4488	Cassette Latch Tool	\$ 224.91
004722TL	TTL4722	Pressure Gauge (ML700)	\$ 47.31
004723TL*	TTL4723	Pressure Gauge Tool (ML700)	\$ 10.00
G9903246	T9903246	Vacuum Message RAET-EXP.	\$ 65.65
004728TL*	TTL4728	Adapter Tool	\$ 10.00
SP29010188	T29010188	AdapterF. Vac. Gauge (ML700)	\$ 50.21
TL4729	TTL4729	Antistatic Wristlet	\$ 46.90
TL4740*			
G9901918	T9901918	Plastic Key	\$ 7.11
TL4804	TTL4804	Gage-Cassette Opener (ML700)	\$583.28
TL4830	TTL4830	Blow Pipe Positioner-Cassette (ML700)	\$120.58
TL5155	TTL5155	Modem Cable (Laser)	\$ 18.73
TL5303	TTL5303	Gloves Neoprene Acid Resist	\$ 15.56
TL5541	TTL5541	Manifold Test Kit-Vacuum	\$162.48
TL5542	TTL5542	Manifold Test Kit-Pressure	\$179.29
TL5543	TTL5543	Coupling-Black	\$53.08
TL5544	TTL5544	Coupling-Red	\$57.91
005335K	T005335K	Standby Timer Kit	\$448.73
N.A.	T0189	Bodine Dual-Speed Motor	\$788.02
1C0030	T1C0030	Main Harness Assembly (M7B)	\$347.51
1C0033	T1C0033	Buzzer Assembly (M7)	\$ 35.41
1C0328	T1C0328	Rack Bar (M6 series)	\$517.80
1C0489	T1C0489	D/F Cross-Over Roller – White Silicone (5000RA)	\$330.56
1C0950	T1C0950	Air Tube (M6 series)	\$ 31.67
1C4315	T1C4315	Pulley (silicone)	\$ 27.78
1C4384	T1C4384	Tube-Yoke Supply (270/3000RA)	\$1,166.55
1C4385	T1C4385	Hose-Yoke (270/3000RA)	\$ 10.00
1C4395*			
1C4398	T1C4398	Tubing	\$ 32.13
1C4564*			
3H4056	T3H4056	Solution Heater (180LP)	\$720.20
1C4593	T1C4593	Replen. Pump (180,270,460,480,3000,5000RA)	\$1,808.15
1C4951	T1C4951	Shaft, Dryer (180LP)	\$ 70.50
1C4952	T1C4952	Drive Shaft, Dryer (270RA)	\$271.53
1C4968	T1C4968	Valve (.75 gpm)	\$114.79
1C4981	T1C4981	Spacer	\$ 5.00
1C7006	T1C7006	Bellows	\$105.40
1C7106	T1C7106	Spring-Extension	\$ 44.54
1C7107	T1C7107	Ext. Spring	\$10.92
1C7237	T1C7237	Shoulder Screw-Level Sensor	\$ 16.99

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
1C7315*			
7C7862	R7C7862	Developer Rack (M35M/M35A-M)	\$3,009.35
1C7316	T1C7316	Clip, Solution-Unpainted (M35 series)	\$159.08
1C7317	T1C7317	Clip-Tank End-Unpainted(M35 series)	\$224.28
1C7318	T1C7318	Sprocket (22-tooth)	\$115.23
1C7321	T1C7321	Mammo Kit (M6B,M6AN,M6AW series)	\$682.03
1C7366	T1C7366	Shaft, Drive 3000RA)	\$118.69
1C7377	R1C7377	Detector Assy. (M6B series-S/N 10,000*)	\$3,137.44
1C7381	T1C7381	Lustran Utility Roller	\$ 36.98
1C7383	T1C7383	Roller, Deflector Assy. (M7 series)	\$ 55.57
1C7384	T1C7384	Lustran "A" Roller	\$124.68
1C7385	T1C7385	Roller	\$ 55.63
1C7386	T1C7386	Lustran "D" Roller	\$ 64.01
1C7387	T1C7387	Lustran "H" Roller	\$ 23.34
1C7389	T1C7389	Lustran "E" Roller	\$ 23.82
1C7390	T1C7390	Lustran Roller	\$ 24.23
1C7391*			
1C7381	T1C7381	Lustran Utility Roller	\$ 36.98
1C7393*			
1C7381	T1C7381	Lustran Utility Roller	\$ 36.98
1C7395	T1C7395	Roller Assy.	
1C7396	T1C7396	Deflector Roller "J" (M7)	\$ 38.24
1C7397	T1C7397	"W" Roller	\$ 31.82
1C7398	T1C7398	Knurled Roller	\$ 36.98
1C7550	T1C7550	Roller	\$223.97
1C7551	T1C7551	Lustran "D" Roller	\$247.91
1C7552	T1C7552	Roller, Detector Assy.	\$70.20
1C7554	T1C7554	"D" Roller (M35 series)	\$ 68.92
1C7555	T1C7555	Lustran Roller	\$ 43.91
1C7556	T1C7556	Lustran Roller w/gear & sprocket (M35 series)	\$ 49.20
1C7557	T1C7557	Knurled Roller (M35 series)	\$ 78.06
1C7558	T1C7558	Lustran "A" Roller (M35 series)	\$ 61.05
1C7630	T1C7630	Sideplate, Adjustable Drive-Side	\$ 33.33
1C7631	T1C7631	Sideplate, Adjustable Non-Drive-Side	\$ 63.63
1C7632	T1C7632	Sideplate	\$ 10.00
1C7633	T1C7633	Rubber "F" Roller (M35 series)	\$137.78
1C7635	T1C7635	Guide Shoe Assy. (M35 series-bottom)	\$107.73
1C7636	T1C7636	Stainless "G" Roller (M35 series)	\$582.31
1C7640	T1C7640	Rubber "B" Roller Assy. (M35 series)	\$319.91
1C7772	T1C7772	Relay	\$327.79
IF5313	T1F5313	Relay	\$228.79
1F5709	T1F5709	Kit-Smooth Guidesheets (M6B/460RA/480RA/5000RA)	\$249.29
N.A.	T110-016I	Internal Standby Harness	\$142.88
N.A.	T110-017E	External Standby Harness	\$127.40
SP126869	T126869	Nut	\$5.00
133589	T133589	Plug	\$47.38
135119	T135119	Receptacle	\$34.71
SPSP136224	T136224	E-Ring (M7B)	\$4.64
141-6957	T141-6957	Standby Control (M6AN/AW-includes solenoid)	\$2,474.62
144097	T144097	Wing Nut	\$ 8.54
SP164385	T164385	Pin	\$7.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
165152	T165152	Metal Jumper	\$5.00
165626	T165626	Screw	\$5.00
165627	T165627	Washer	\$5.00
SP165678	T165678	Standby Screw	\$3.87
165938	T165938	Clamp <i>(with nut)</i>	\$11.49
166323	T166323	Pin	\$5.00
166326	T166326	E-Ring	\$5.00
5B6850*			
166331	T166331	Nut	\$5.00
166856	T166856	Screw	\$5.00
167564	T167564	Cover	\$127.89
168424	T168424	Washer-Lock	\$ 5.00
169683	T169683	Pin	\$.90
N.A.	T17494	Capacitor <i>(for new-style M6 Main Dr. Motor)</i>	\$98.41
175381	T175381	Pin	\$5.00
SPK175760	T175760	Retaining Ring	\$.52
SPK176250	T176250	Pin	\$3.09
177006	T177006	Potentiometer	\$30.27
178107	T178107	Screw	\$5.00
179720	T179720	Connector	\$ 43.91
180281	T180281	Connector-Tee	\$5.00
180557	T180557	Tee Connector	\$ 63.63
180977	T180977	Washer	\$5.00
181587	T181587	Relief	\$.76
181637*			
609952	T609952	Flow Control Valve <i>(1.5 gpm)</i>	\$142.19
SP182590	T182590	Tie-Cable	\$1.62
SPKC182915	T182915	Tie Cable <i>(8")</i>	\$1.03
183153	T183153	Heater <i>(M6AN/AW series)</i>	\$254.52
183302	T183302	Screw	\$1.60
183322	T183322	Screw	\$5.00
183624	T183624	Screw	\$5.00
183851	T183851	O-Ring	\$ 5.00
183882	T183882	Fitting, Hose	\$5.10
184815*			
914443	T914443	Circuit Breaker-30amp <i>(M6 series)</i>	\$205.51
184816	T184816	Block-Terminal	\$54.05
184823	T184823	Nipple – Close	\$10.00
185101	T185101	Nipple – Close <i>(3/8" NPT)</i>	\$12.76
185480	T185480	Lamp	\$34.34
185896	T185896	Thumb Screw <i>(M6 series)</i>	\$13.45
186153	T186153	Switch	\$ 69.74
186441	T186441	Tee	\$4.05
187248*			
914444	T914444	Circuit Breaker, 20 amp <i>(M7B)</i>	\$185.79
187283	T187283	Lockwasher	\$ 5.00
187296	T187296	Terminal Strip	\$127.89
187302*			
183882	T183882	Adapter	\$5.10
187632	T187632	Adapter-Hose	\$5.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
187774	T187774	Screw	\$5.00
188681*			
246842	T246842	K6 Relay (3-sec.)	\$590.88
188967	T188967	Link, Master	\$5.00
189734	T189734	Receptacle	\$6.83
189769	T189769	Small Standby Relay-K3	\$84.29
189770	T189770	Relay	\$167.96
189771*			
189772*			
7C8478	T7C8478	Water Solenoid (M6 series)	\$3,134.32
189773	T189773	Bushing	\$5.00
190481	T190481	Bearing (small)	\$10.00
190596	T190596	Screw	\$5.00
191097	T191097	Screw	\$5.00
191098	T191098	Washer, Lock #4, Stainless	\$5.00
192019	T192019	Bearing, Large	\$53.80
192203	T192203	Screw (M7B feed tray stand)	\$2.15
194116	T194116	Timing Belt	\$24.16
2E4319	T2E4319	Guard, Splash: Stain	\$205.57
2E4648	T2E4648	Relay (K1, K2, K3)	\$87.51
240002	T240002	Tube, Dev/Replen.	\$20.68
240003	T240003	Tube, Fix/Replen.	\$20.68
240004	T240004	Formed Tube	\$29.35
240005	T240005	Formed Tube	\$27.91
240058	T240058	Tie-Red	\$5.00
240059	T240059	Tie-Blue	\$5.00
240060	T240060	Bearing—Plain	\$9.90
240061	T240061	Bearing—Flat	\$13.65
240065	T240065	Sleeve Spacer	\$10.00
240066	T240066	Roller—Crossover (Drive)	\$352.17
240067	T240067	Roller-Crossover (Driven)	\$371.20
240105	T240105	Tubing, Formed 72 in. (M6B series)	\$21.23
240106	T240106	Tubing, Formed 69 in. (M6B series)	\$15.63
240109	T240109	Detector Switch	\$26.31
240110	T240110	Det. Switch Cable (M6B series)	\$95.45
240111	T240111	Formed Tubing 64 in. (M6B series)	\$30.64
240113	T240113	Cable	\$44.54
240116	T240116	ED Cable (M6B series)	\$77.30
240117	T240117	ED-Jumper XMFR	\$10.00
240225	T240225	Blower Assy. (2000A series)	\$211.05
240227	T240227	Pump, Replen.-115V (2000A series)	\$2,041.46
240235	T240235	Bin-Receiving	\$961.52
240416	T240416	Insulation, Black (per sheet)	\$46.54
240514	T240514	Film Detector Switch (M35A/AM series)	\$81.46
240566	T240566	Nameplate(2000A series))	\$194.23
240627	T240627	Terminal Strip	\$197.69
240628	T240628	Tubing, Formed (M6B series)	\$12.56
240634	T240634	Tubing, Formed 15.5" (M6B series)	\$9.49
240680	R240680	Developer Rack Assembly	\$5,674.96
N.A.	S240680	New Style Dev. Rack Assy. (for Mammo)	\$5,896.29
240681	R240681	Fixer Rack Assembly (M6AN/AW series)	\$5,396.29

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
240682	R240682	Wash Rack Assembly (M6AW & M20 series)	\$5,719.96
240683*			
1C7377	R1C7377	Detector Assembly 3-roller (M6AN/AW, M20, M6B)	\$3,137.44
240683	R240683	Detector Assembly-2 roller (M6AN/AW, M20, M6B)	\$1,485.00
240684	R240684	Dev/Fix X-Over Assy. (M6B/480RA series – Silicone)	\$972.80
N.A.	S240684	Dev/Fix X-Over Assy. (M6B/480RA series-mammo)	\$1,736.37
240685*			
779155	R779155	Fix/Wash X-Over Assy. (M6AN/AW series)	\$1,004.73
240686	R240686	Squeegee Assembly (M6AN/AW series)	\$3,162.27
240726	T240726	Protector & Shelf Assembly	\$193.85
240735	T240735	Sprocket, 20-tooth	\$98.09
240747	T240747	Developer/Fix Cross-Over Roller (M35/M35M)	\$352.42
240819	T240819	Drain Valve	\$66.31
240868	T240868	Spacer, Electrical Box (M6B series)	\$44.98
240888	T240888	Valve, Replen.	\$92.02
240894	T240894	Bracket, Thermowell (M6 series)	\$143.77
N.A.	R240901-AN	Reconditioned* Stainless Tank (M6A/N series)	\$4,001.01
N.A.	R240901-AW	Reconditioned* Stainless Tank (M6AW series)	\$4,001.01
240912	T240912	Cable Assy. (M6B series)	\$163.28
240916	T240916	Cable Assy., Dryer	\$173.75
246531*			
240235	T240235	Film Receiving Bin	\$961.52
246554	T246554	Adapter Plate (use with T914441 circuit breaker)	\$507.15
246575	T246575	Roller (M6I)	\$614.29
246576	T246576	Roller (M6I)	\$447.49
246579	R246579	Dev/Fix X-Over Assembly (M6-I series)	\$1,888.03
246580	T246580	Spring (Min-R)	\$11.12
246581	T246581	Gear	\$11.47
246588	T246588	Sprocket, 12-tooth	\$116.42
246590	T246590	Sprocket, 12-tooth (with 2 setscrews)	70.43
246612*			
1C0328	T1C0328	Rack Bar (M6 series)	\$517.80
246747	T246747	Blower Motor Guard	\$582.24
246800	T246800	Clamp (3/8")	\$5.00
246801	T246801	Clamp (1/2")	\$5.00
246802*			
246801	T246801	Clamp (1/2")	\$5.00
246803	T246803	Clamp (3/4")	\$5.00
246804	T246804	Clamp	\$5.00
246806	T246806	Clamp, Exhaust Hose	\$5.00
246808	T246808	Clamp (3/4")	\$10.00
246840	T246840	Electrical Box (M35A series)	\$2,949.45
246842	T246842	K6 Relay (M6B series)	\$590.88
246868	T246868	Relay Kit (M7B, M35, M35A)	\$415.49
246879	T246879	Sideplate (entrance)	\$8.94
246882	T246882	Upper Detector Roller (M35 series)	\$287.66
246883	T246883	Lower Detector Roller (M35 series)	\$205.25
246890	T246890	Tube	\$46.75
246911	T246911	Drip Guard	\$90.47
246915	T246915	Roller Assy.-squeegee (Silicone-Wash Rack)	\$371.07
246916	T246916	Roller Assy.-squeegee (Silicone-Wash Rack)	\$386.19

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
246917	T246917	Silicone "G" Roller (M35AM/M series)	\$423.49
246918	T246918	Silicone "F" Roller (M35AM/M series)	\$365.78
246963	T246963	Harness Assy.	\$174.32
246974	T246974	Tubing, Formed (M6B series)	\$23.34
246979	T246979	Roller, Phenolic	\$103.64
246982	T246982	Shaft-Driven (270RA wash)	\$333.14
261180*			
7B8229	T7B8229	Captive Screw	\$7.85
261275	R261275	Developer Rack Assembly (M35 series)	\$4,110.88
261276	R261276	Fixer/Wash Cross-Over (M35 series)	\$1,606.20
261277	R261277	Wash Rack Assembly (M35 series)	\$4,883.17
261279*			
1C7384	T1C7384	Lustran "A" Roller	\$124.68
261284*			
651125	T651125	Belt, Dryer Transport (M6 series)	\$80.01
261292	T261292	Blower Wheel (M6B series)	\$91.85
261376	T261376	Valve & Tee Assembly	\$73.71
261385	T261385	Tubing, Flared (M6B series)	\$9.69
261386	T261386	Orifice (.187 diameter)	\$62.69
261389	T261389	Screw (M6B tank)	\$5.00
261590*			
261650	T261650	Flow Control Valve Kit	\$382.91
261632	T261632	Back Panel (M7B)	\$158.95
261650	T261650	Flow Control Valve Kit	\$382.91
261664	T261664	Plug (M6 tank)	\$37.67
261784*			
3H9103	T3H9103	Support Assy.-Feed End	\$31.82
261885	T261885	M6B PC Board (S/N 10,000 & above)	\$2,167.94
261895	R261895	Wash Rack Assembly (460/480RA series)	\$1,970.51
261896	R261896	Fixer Rack Assembly (460/480RA series)	\$4,524.86
261897	R261897	Fix/Wash X-Over Assembly (460/480RA & M6 series)	\$5,671.17
261898	T261898	Cup, for leveling screw	\$12.72
264405	T264405	Fuse Block (M6B series)	\$46.78
264407	T264407	Fuse	\$59.08
264425	T264425	Exhaust Hose	\$31.69
264426	T264426	Plate, Motor Mount	\$521.77
264445	T264445	Top Panel Assembly (M35 series)	\$411.71
264482	T264482	Dryer Roller (M7 series)	\$113.27
264483	T264483	Sprocket, 22-tooth	\$66.25
264484	T264484	Cooling Solenoid (M6B series)	\$130.73
264490	T264490	Cable Assembly	\$353.75
264503	T264503	Exhaust Fan Kit, 110VAC-60Gz	\$2,249.30
264519	T264519	Duct Assy.-Vacuum Adjust.	\$940.23
264526	T264526	Splash Guard (M6 series)	\$602.09
264528*			
264539	T264539	Dryer Plenum (M6 series)	\$545.96
264529*			
7C6650	T7C6650	Recirc. Pump Guard	\$74.28
264530	T264530	Drip Guard (M6 series)	\$183.14

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
264531*			
246747	T246747	Blower Motor Guard	\$582.24
264533	T264533	Guard, Chain	\$56.76
264535	T264535	Filter Assy., Drip Pan	\$134.38
264536	T264536	Splash Guard (M6 series)	\$197.57
264537	T264537	Guard	\$51.42
264539	T264539	Dryer Plenum (M6 series)	\$545.96
264545	T264545	Lip Assy., Drip Cover (M6 series)	\$212.84
264762	T264762	Lamps, DS1/DS2	\$34.67
264832	T264832	Bracket, Thermowell	\$54.94
264834	T264834	Cable Assy.-Solenoid (M6B -[10,000 ^{+D}])	\$58.03
264842	T264842	Front Panel-no bin or bezel included! (M6 series)	\$1,487.63
264843	T264843	Top Cover Assy. (M6 series)	\$1,888.14
264844	T264844	Side Panel (M6 series)	\$1,532.13
286900	T286900	Sideplate, Crossover	\$10.00
286902*			
911809	T911809	Sideplate, Fixer (DS-M7)	\$25.32
286903*			
286904	T286904	Sideplate, Developer (NDS-M6 series)	\$29.14
286907	T286907	Sideplate, Developer (DS-M6 series)	\$33.31
286910	T286910	Sideplate, Wash (NDS-M6 series)	\$25.46
286912	T286912	Sideplate, Wash (DS-M6 series)	\$24.91
286914	T286914	Sideplate, Fixer (NDS-M6 series)	\$27.30
286917	T286917	Sideplate, Fixer (DS-M6 series)	\$26.96
286920*			
911845	T911845	Sideplate, Fixer (NDS-M7 series)	\$25.25
286921*			
286922	T286922	Sideplate, Squeegee (DS/NDS)	\$12.97
286924	T286924	Sideplate, Wash (DS-M7 series)	\$26.89
286926	T286926	Sideplate, Wash (NDS-M7 series)	\$27.71
286928	T286928	Sideplate, Wash (NDS-M6 series)	\$44.16
286930	T286930	Sideplate, Wash (DS-M6 series)	\$68.87
286933	T286933	Sideplate, Developer (DS-M7 series)	\$43.31
286934	T286934	Sideplate, Developer (NDS-M7 series)	\$94.12
286948	T286948	"E" Roller Assy.	\$175.64
286949	T286949	"C" Roller	\$257.04
287002*			
7C6660	T7C6660	Tie Rod (M35 series)	\$25.59
287003	T287003	Shaft (M35 series)	\$33.58
287004*			
5B6077	T5B6077	"B" Roller Assembly (M35 series)	\$233.98
287005	T287005	Rubber "A" Roller (M35 series)	\$201.03
287009*			
1C7635	T1C7635	Bottom Turnaround Guide Shoe (M35 series)	\$107.73
287013*			
246917	T246917	Silicone "G" Roller (M35 series)	\$423.49
287015*			
246918	T246918	Silicone "F" Roller (M35 series)	\$365.78
287018	T287018	Guide Rod Assy. (M35 series)	\$170.54
287020*			
1C7554	T1C7554	"D" Roller (M35 series)	\$68.92

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
287021*			
1C7555	T1C7555	Lustran Roller	\$43.91
287025*			
261277	R261277	Wash Rack Assy. (M35 series)	\$4,883.17
287027	T287027	Drive Shaft (M35 series)	\$55.25
287028	T287028	"G" Roller (M35 series)	\$120.52
287030	T287030	"F" Roller (M35 series)	\$581.62
287032	T287032	Rewet Roller (M35 series)	\$169.91
287035	T287035	Tie Rod, Detector Assy. (M35 series)	\$38.43
287036*			
3E0999	T3E0999	Shaft, Tie Rod Det. Assy. (M35 series)	\$34.52
287037	T287037	Shaft-Entrance (M35 series)	\$71.25
287039	T287039	Shaft (M35 series)	\$41.27
287041	T287041	Textured Entrance Guides hoe (M35 series)	\$58.84
287047	T287047	Shaft-Drive (M35 series)	\$81.77
287048	T287048	Rubber Cross-Over Roller (M35 series)	\$123.17
287052*			
748147	T748147	Cross-Over Roller (M35 series)	\$121.91
287053		Heat Exchanger (M35 series)	\$134.06
287059*			
1C7555	T1C7555	Lustran Roller	\$43.91
287066*			
1C7556	T1C7556	Lustran Roller-w/gear & sprocket (M35 series)	\$49.20
287067	T287067	Roller, Dryer (Phenolic) (M35 series)	\$102.25
287078	T287078	Roller, Rubber (M35 series)	\$109.18
287079	T287079	Shaft, Upper Roller (M35 series)	\$118.82
287081	T287081	Roller, Dryer (Phenolic) (M35 series)	\$182.26
287085	T287085	Shaft, Lower Roller (M35 series)	\$20.41
287088	T287088	Shaft, Dryer Drive (M35 series)	\$28.12
287656	T287656	Roller, Dryer-w/o bearing (M35 series)	\$161.34
287676	T287676	Tie Rod (M35 series)	\$20.95
287677*			
1C7316	T1C7316	Clip-Solution -Unpainted (M35 series)	\$159.08
287679*			
1C7317	T1C7317	End Clip -Unpainted (M35 series)	\$224.28
287680	T287680	Cover-Evaporation	\$34.34
287683	T287683	Seal	\$5.00
287685	T287685	Feed Tray Overlay Protector (M35 series)	\$105.65
287689	T287689	Tie Rod (M35 series)	\$34.34
287730	T287730	Water Bracket Assy.	\$43.22
287733	T287733	Tubing (M35 series)	\$32.76
287735	T287735	Tubing (M35 series)	\$18.91
287736	T287736	Tubing (M35 series)	\$39.25
287739*			
914568	T914568	Tube-Filter/Thermowell	\$32.76
287740	T287740	Tubing (M35 series)	\$30.23
288210	T288210	Cable, J5/P5	\$116.36
288211	T288211	Cord-Dryer Heater (M6B series)	\$258.49
288215	T288215	Dryer Heater Assy. (M6B series)	\$852.37
288233	T288233	Line Filter (M6B series)	\$99.16
289783	T289783	Pulley/Bearing Assy.	\$56.67

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
3B7301	T3B7301	Thrust Bearing	\$12.22
3E0957	T3E0957	Shaft, Main Drive (Min-R)	\$68.48
3E0958	T3E0958	Bearing, Main Drive (Min-R)	\$16.18
3E0959	T3E0959	Spacer, Rack Support, Drive Side (Min-R)	\$32.51
3E0961	T3E0961	Spacer, Rack Support, Non-Drive Side (Min-R)	\$6.89
3E0962	T3E0962	Tie Rod, Developer (Min-R)	\$30.03
3E0967	T3E0967	Bearing, Dev. Rack (Min-R)	\$105.27
3E0976	T3E0976	Roller, Entrance, Dev. Rack (Min-R)	\$2.76
3E0979	T3E0979	Snap-On Spacer (Min-R)	\$15.83
3E0982	T3E0982	Shaft, Master Roller (Min-R)	\$215.84
3E0984	T3E0984	Shaft, Entrance (Min-R)	\$22.93
3E0989	T3E0989	Washer, Rewet (Min-R)	\$112.33
3E0990	T3E0990	Bearing, Rewet (Min-R)	\$10.17
3E0995	T3E0995	Spring-Exit Pinch (Min-R)	\$78.56
3E0996	T3E0996	Spring-Turnaround-V (Min-R)	\$41.33
3E0997	T3E0997	Spring-Turnaround-H (Min-R)	\$5.00
3E0999	T3E0999	Tie Rod, Detector (Min-R)	\$34.52
3E1938	T3E1938	Roller, Drive, D/F X-Over (Min-R)	\$364.20
3E1939	T3E1939	Roller, Driven, D/F X-Over (Min-R)	\$296.60
3E1943	T3E1943	Dev/Fix X-Over Assy. (Min-R)	\$2,505.02
3E1944	T3E1944	Min-R Developer Rack	\$5,815.74
3E1945	T3E1945	Detector Assy. (Min-R)	\$2,138.66
3H0648*			
3H3936*			
4E8210	T4E8210	Quad DC Power Supply	\$575.57
3H4056	T3H4056	Solution Heater (180LP)	\$720.20
3H4125	T3H4125	Sideplate, Entrance Detector (Min-R)	\$7.99
3H4146	T3H4146	Sideplate, D/F X-Over (Min-R)	\$20.13
3H4475	T3H4475	Main Drive Motor (M6A series)	\$1,098.29
3H8726	R3H8726	Fix/Wash X-Over (5000RA)	\$3,566.93
3H8738	T3H8738	Chain - Drive (32 links)	\$122.03
3H9103	T3H9103	Support Assy.-Feed End	\$31.82
3H9104	T3H9104	Support Casting-Dryer End	\$575.63
333428	T333428	Tubing (3/4" I.D.)	\$10.00
333432	T333432	Tubing, Formed (M6B series)	\$28.80
333437	T333437	Tubing(M35 series)	\$28.53
333438	T333438	Flared Tube (M35 series)	\$31.82
333439	T333439	Tubing, Flared(M35 series)	\$12.29
333440	T333440	Tubing, Flared(M35 series)	\$63.63
333441	T333441	Tubing, Formed(M35 series)	\$10.72
333452	T333452	Tubing, Formed (M6B series)	\$9.56
333464	T333464	Tubing, Formed(M35 series)	\$15.15
333465	T333465	Tube (M6B series)	\$27.91
333519	T333519	Pump, Recirc.- 220V (M35/M series)	\$1,192.77
333543*			
532054	T532054	Bearing (M35 dryer)	\$7.78
333599	R333599	Reconditioned Input Roller (autofeeder)	\$532.41
333600	R333600	Reconditioned Output Roller (autofeeder)	\$443.96
333658	T333658	Main Drive Motor (M6B -S/N above 10,000)	\$1,093.03
333659	T333659	Pump, Replen.-240V (M6B -S/N above 10,000)	\$2,303.94
333789	R333789	Turnaround Assy. (w/T5B6076 "B" roller)	\$947.16

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
333794*			
5B6163	R5B6163	Detector Entrance Assy. (M35 series)	\$2,164.76
333795	R333795	Developer/Fixer Cross-Over-(Silicone) (M35A)	\$980.77
N.A.	S333795	Rubber Developer/Fixer X-Over (M35A series)	\$1,489.89
333961	T333961	Nozzle Clip, PHII	\$10.00
333962	T333962	Tubing, Formed (M6B series) PHII	\$18.22
333966	T333966	Tubing, Formed	\$15.29
333973	T333973	Tubing, Formed	\$18.91
333975	T333975	Access Cover (M35 series)	\$98.66
333978	T333978	Front Panel (M35 series)	\$210.48
333981	T333981	Side Panel (M35 series)	\$249.10
333994	T333994	Exit Panel (M35 series)	\$289.30
335534	R335534	Fixer Rack (480RA)	\$5,244.79
335656*			
9B8088	T9B8088	Squeegee Kit/Driven Roller	\$467.71
335657*			
9B8087	T9B8087	Squeegee Kit/Drive Roller	\$334.22
335663	T335663	Shaft, Sprocket Dryer (270RA Integrated)	\$106.91
335975	R335975	Developer/Fixer Cross-Over (M7B)	\$1,159.04
N.A.	S335975	Silicone Developer/Fixer X-Over (M7B)	\$2,100.00
335976	R335976	Fixer/Wash Cross-Over (M7B)	\$1,837.43
335977	R335977	Turnaround Assy. (w/T536884 rubber "B" roller)	\$2,927.29
4B0985	T4B0985	Line Filter	\$43.16
4E8210	T4E8210	Quad DC Power Supply	\$575.57
418700	T418700	Impeller Housing	\$114.53
418701	T418701	Housing & Impeller Assembly	\$292.70
418702	T418702	Pump, Recirc.-230V (M6AN/AW)	\$1,908.46
418703*			
333519	T333519	Pump, Recirc.- 220V (M35/M series)	\$1,192.77
418709	T418709	Pump, Recirc.-115V (M7B)	\$797.35
418724*			
650852	R650852	<i>Reconditioned</i> Developer Thermostat	\$2,000.17
N.A.	R650852	<i>Reconditioned</i> Developer Thermostat (<i>without thermistor!</i>)	\$1,600.50
418726	T418726	Variable Resistor-POT	\$119.57
418729*			
185480	T185480	ED-Lamp	\$34.34
418730	T418730	Thermistor, 50K (M6 series)	\$1,726.90
418734	T418734	Solid State Relay	\$372.58
418746	T418746	Main Harness Assembly (M6 series)	\$1,688.43
N.A.	T418746-w/o	Harness Assembly (M6 series)-(without breakers)	\$639.00
418775	T418775	Washer-Thrust	\$5.00
419788*			
1C7557	T1C7557	Knurled Roller (M35 series)	\$78.06
419779	T419779	Flooded Replenishment Kit (M35 series)	\$2,714.37
419780*			
7C6980	T7C6980	Replenishment Timer (M35 series)	\$439.05
420224	T420224	Connector (M6B series)	\$7.31
422206	T422206	Jack, Connector	\$10.00
431908	T431908	Pin	\$16.93
452128*			
452044	T452044	Key, long	\$7.51

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
452101	T452101	Nut	\$5.00
452111	T452111	Ladder Chain	\$95.05
452673	T452673	Filter	\$23.17
452708	T452708	Tubing Bushing	\$5.00
452967	T452967	Bell	\$539.85
N.A.	S452967	M6A <i>Electronic Buzzer</i> (the bell alternative)	\$539.85
452990	T452990	Tubing (3/8") (sold per ft.)	\$10.00
454260	T454260	Master Link	\$5.85
454285	T454285	Drive Sprocket	\$10.00
454354	T454354	Key, short	\$5.00
454448	R454448	<i>Reconditioned</i> Main Drive Motor (Ind. B)	\$2,530.63
454454	T454454	Sprocket, Motor	\$48.39
454543	T454543	O-Ring	\$5.00
454648	T454648	Grommet (small)	\$.83
454749	T454749	Gasket(per14-ftlength)	\$22.38
454984	T454984	Grommet (medium)	\$1.05
456162	T456162	Shaft, Dryer (M6 series)	\$33.39
456560	T456560	Nipple(M6Bseries)	\$10.00
457253	T457253	Processor Panel Gasket (per foot)	\$2.63
457672	T457672	Bearing, Blower Upper	\$101.04
458161*			
650705	T650705	Bearing	\$11.06
458164	T458164	Shaft	\$79.19
N.A.	TL458164	Shaft(M6series,mammocross-over)	\$109.05
458170	T458170	Tie Rod (M6 series)	\$41.33
458174	T458174	Sleeve Spacer	\$10.00
458175	T458175	Shaft (M6 series)	\$53.49
458182*			
550781	T550781	Gear, X-Over	\$10.03
458189*			
748955	T748955	Worm Gear, L.H. (□White□)	\$10.00
458190*			
286922	T286922	Sideplate (squeegee)	\$12.97
458191*			
748956	T748956	Worm Gear, R.H. (□Black□)	\$10.00
458205	T458205	Spacer	\$5.00
458207	T458207	Shaft	\$45.80
458209	T458209	Spacer (for old-style sideplate)	\$10.65
458210*			
559787	T559787	Shaft	\$51.91
458219*			
7C6662	T7C6662	Tie Rod (M6/M7 series)	\$24.50
458235	T458235	Shaft, Entrance (M6 series)	\$49.27
458236	T458236	Small Sleeve Spacer	\$19.59
458240	T458240	Chain, Wash (M6 series)	\$48.07
458350	T458350	Feed Stand & Guide (M6 series)	\$98.66
458382	T458382	Gear, Spur	\$27.57
458383*			
650663	T650663	Drive Gear	\$9.49
458384	T458384	Gear	\$10.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
458386	T458386	Screw	\$47.06
458387*			
1C4981	T1C4981	Spacer	\$5.00
458389	T458389	Bearing, dryer squeegee	\$22.32
458390	T458390	Shaft, Dryer Squeegee (M6series)	\$31.46
458391	T458391	Shaft (M6 squeegee, "J" Roller)	\$33.24
458393	T458393	Shaft (M6 squeegee, "L" Roller)	\$76.92
458394*			
7C6661	T7C6661	Tie Rod (M6 series-entrance)	\$32.46
458423	T458423	Squeegee "J" Roller	\$140.30
458505	T458505	Bracket	\$37.45
458508	T458508	Pulley-Belt (M6AN/AW series)	\$71.88
458610	T458610	Bezel (M6 series)	\$107.92
458709	T458709	Support Block	\$5.00
458850	T458850	Exhaust Hose (2-ply; 2-inch)	\$20.27
458901*			
1C4315	T1C4315	Pulley (silicone)	\$27.78
458946	T458946	Pulley (M6A series)	\$41.33
458949	T458949	"N" Squeegee Roller (M6 series)	\$144.14
458956	T458956	Belt	\$100.43
458977	T458977	Button, Plug	\$20.61
458979	T458979	Check Valve	\$381.78
458980*			
581200	T581200	Drain Valve	\$92.67
459093	T459093	Ball Float	\$87.21
459229	T459229	Leveling Screw	\$5.00
459248	T459248	Check Valve (M6 series)	\$32.28
459263	T459263	Chain, Developer (M6 series)	\$56.45
459333	T459333	Chain Segment	\$11.94
459343	T459343	Washer	\$5.00
459344*			
1C7106	T1C7106	Spring-Extension	\$44.54
459346	T459346	Spacer	\$7.11
459347	T459347	Collar	\$19.79
459354	T459354	Washer, Spl.	\$.75
459355	T459355	Spacer	\$5.00
459378	T459378	Tie Rod with handle (M6 series)	\$73.84
459379	T459379	"C" Roller (M6 series)	\$133.94
459473	T459473	Knob for Drain Valve	\$18.06
459548*			
261898	T261898	Plate, Floor	\$12.72
459576	T459576	Orifice, Fixer Recirc.	\$23.96
459600	T459600	Tubing	\$15.49
459609	T459609	Tube, Drain	\$5.00
459700	T459700	Washer	\$19.72
459943	T459943	Ident. Washer-Fixer Rack	\$5.00
459944	T459944	Ident. Washer-Fixer Rack	\$5.00
459946	T459946	Washer, Fixer Rack	\$5.00
459970	T459970	Connector	\$5.00
459990	T459990	Nut-Drain Valve	\$26.07
459991	T459991	O-Ring	\$8.67

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
459992	T459992	Stem & Disk-Drain Valve	\$49.14
459993	T459993	Seal Ring for Drain Valve 4	\$10.37
459994	T459994	Drain Valve Body	\$44.73
459999	T459999	Cover, Bell (M6 series)	\$41.77
460422*			
333789	R333789	Turnaround Assy. (w/T5B6076 "B" roller)	\$947.16
460674*			
418746	T418746	Main Wiring Harness (M6 series)	\$1,688.43
460784*			
474599	T474599	Rewet Roller (M6/M7 series)	\$218.17
460937*			
462048	T462048	Dryer OT (M6A series)	\$13.49
462139	T462139	Washer, "E" Roller	\$10.00
462452*			
651125	T651125	Belt, Dryer Transport (M6 series)	\$80.01
462602	T462602	Tubing (18")	\$15.77
462688	T462688	Chain, Fixer (M6 series)	\$60.04
462719	T462719	Tube, Formed (M7 series)	\$10.00
462885	T462885	Tee	\$10.00
463004	T463004	Spacer-BRG	\$5.00
463009	T463009	Nipple-Drain Valve	\$30.99
463100	T463100	Washer	\$5.00
463363	T463363	O-Ring, Impeller (large)	\$20.68
463950	T463950	O-Ring	\$5.00
464096*			
286900	T286900	Sideplate, Crossover	\$10.00
464101	T464101	Tube	\$10.51
464936	T464936	Elbow 90°	\$10.51
464937*			
180557	T180557	Tee Connector	\$63.63
465361*			
465392	T465392	Chain, Main	\$28.13
466556	T466556	Gudgeon Gear	\$5.00
466640	T466640	Sprocket, 17-tooth/2 set screws	\$389.78
466641	T466641	Support	\$61.99
467015	T467015	"F" Roller, Stainless Steel (M6/M7 series)	\$580.80
467016	T467016	"K" Roller	\$100.86
467017	T467017	"E," "P," "R" Roller	\$124.99
467555*			
1C7387	T1C7387	Lustran "H" Roller	\$23.34
467562*			
1C7389	T1C7389	Lustran "E" Roller	\$23.82
467563	T467563	Gudgeon, Plain Idler	\$13.45
467564	T467564	Gudgeon Gear	\$10.00
467565*			
1C7381	T1C7381	Lustran Utility Roller	\$36.98
467579	T467579	Arm-Take-Up	\$359.10
467843	T467843	Spring	\$10.00
467844*			
474645	T474645	Spring	\$10.00
467867	R467867	Reconditioned Drain Header Assembly	\$1,720.66

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
467919	T482375	Motor Bracket/Angle Plate	\$243.24
467988*			
748437	U748437	Dryer Pan	\$283.25
467990	T467990	Orifice (.250ID)	\$128.33
468062	T468062	Tie Rod with handle (M6 series-entrance & squeegee)	\$286.34
468218	T468218	Receptacle	\$90.91
468220	T468220	Front Plate (w/o temperature cut-out)	\$349.17
468227	T468227	Blower Panel Cover Plate (w/o latch)	\$232.03
468312	T468312	Terminal Strip	\$127.89
468326*			
7C6650	T7C6650	Recirc. Pump Guard	\$74.28
468332	T468332	Pulley-Blower	\$154.74
468354*			
550774	T550774	Gear, Rack	\$7.03
468364	T468364	Mechanical Switch	\$30.32
468521	T468521	Developer Heater	\$506.21
468558	T468558	Formed Tube, PHII	\$10.00
468559	T468559	Tube-Formed Wash	\$8.05
468561	T468561	Tube-Formed Dev.	\$18.97
468564	T468564	Tubing (16.5 in.)	\$7.23
468570	T468570	Formed Tubing	\$9.96
468571	T468571	Formed Tubing	\$16.99
468572	T468572	Formed Tubing	\$63.63
468573	T468573	Formed Tubing	\$10.85
468574	T468574	Cord Assembly	\$150.73
468684	T468684	ED-Jumper XMFR (Red)	\$31.82
468685*			
240117	T240117	ED-Jumper XMFR	\$10.00
468686	T468686	ED-Jumper XMFR (Yellow)	\$18.43
468698	T468698	Latch (M6 series)	\$12.42
468767	T468767	Belt Take-Up Assy.	\$501.10
468967	T468967	Shaft ("B" Roller)	\$374.28
469002	T469002	Screw (each)	\$5.00
469069	T469069	Bearing, Main Drive	\$16.18
469070	T469070	Pulley Assy. (M6 series)	\$130.79
469079	T469079	Bearing (M6 dryer)	\$11.06
469080	T469080	Pulley (M6 series)	\$103.57
469648*			
650711	T650711	Bearing	\$11.12
469349	T469349	22-Link Chain (includes master link)	\$12.01
469401	T469401	Pad	\$20.68
469403	T469403	Pad	\$5.00
469404	T469404	Pad	\$5.00
469409	T469409	Tubing, Formed, Fix.	\$8.53
469411	T469411	Tube, Dev.	\$19.66
469692	T469692	Window Assy. (M6B series)	\$101.49
469702	T469702	Ring-Retaining (each)	\$5.00
470085	T470085	Roller, Rubber Crossover (driven)	\$121.91
470086	T470086	Rubber "G" Roller	\$126.57
470106	T470106	Terminal Strip	\$64.89
470115	T470115	Metering Valve (M6 series)	\$336.67

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
470119	T470119	Handle, Gray	\$6.08
470153	T470153	Handle (M6 series)	\$124.84
470203	T470203	Bearing Assembly	\$258.80
470205	T470205	Blower Shaft (M6 series)	\$534.49
470213	T470213	Bonnet (M6 series)	\$34.65
470214	T470214	Stem (M6 series)	\$75.10
470215	T470215	Gasket (M6 series)	\$8.80
470216	T470216	Body	\$37.55
470465	T470465	"T" Roller (Rubber)	\$164.79
470466	T470466	"T" Roller	\$154.42
470467	T470467	Roller, Rubber "A"	\$152.33
470468	T470468	Rewet Roller	\$608.86
470478	T470478	Harness, Transformer (M6 series)	\$221.95
470519	T470519	Nut, Jam (M6 series)	\$16.58
470534	T470534	Shoulder Screw	\$133.37
470535	T470535	Shoulder Screw	\$105.84
470536	T470536	Shoulder Screw	\$200.84
470581*			
477029	T477029	Clip	\$5.00
470767	T470767	Gear (for T470085 roller)	\$212.06
471335	T471335	Base Plate, Standby Assy. (M6 series)	\$179.47
471347	T471347	Cover, Standby Assy. (M6 series)	\$244.25
471393	T471393	Tube	\$10.00
471397*			
1C0950	T1C0950	Air Tube (M6 series)	\$31.67
472139	T472139	Tubing	\$9.08
472265*			
9B8960	T9B8960	Locking Plate	\$7.03
472223*			
3H9104	T3H9104	Support Casting-Dryer End	\$575.63
472267*			
3H9103	T3H9103	Support Assy-Feed End	\$31.82
472268	T472268	Clip, Cross-Over	\$13.24
472315		Shaft	\$47.12
472337*			
650707	T650707	Bearing	\$10.00
472446	T472446	Drive Shaft	\$188.94
472477	T472477	Flow Control Valve (.25 gpm)	\$155.86
472592*			
5B6527	T5B6527	Bearing, Detector Assembly	\$12.29
472605	T472605	Pump Housing, LH	\$50.90
472606	T472606	Pump Housing, RH	\$50.90
472610*			
418702	T418702	Recirculation Pump-230V (M6A series)	\$1,908.46
472755	T472755	Nut, Tie-Rod (Nylon)	\$10.00
472756	T472756	Nut, Tie-Rod (Nylon)	\$9.21
472761	T472761	Dryer Thermostat(M6series)	\$337.68
472764	T472764	Dryer Knob(M6series)	\$98.66
472811	T472811	Plate, Locking	\$129.28
472914	T472914	Thermostat Cover	\$137.09
472989	T472989	Collar	\$6.83

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
472992*			
472761	T472761	Dryer Thermostat(M6 series)	\$337.68
472994	T472994	Plate	\$139.99
472995	T472995	Plate	\$76.36
473117	T473117	Strap-Filter Assy.	\$127.89
473122*			
264535	T264535	Drip Pan, Filter Assy.	\$134.38
473125	T473125	Formed Tubing	\$22.93
473146*			
5B6076	T5B6076	"B" Roller	\$252.44
473221	T473221	Keeper Plate	\$127.70
473246	T473246	Plastic "T" Connector	\$6.50
473434	T473434	Spacer	\$12.35
473538	T473538	Spacer, Rack	\$10.00
473539	T473539	Spacer	\$27.57
473548*			
286917	T286917	Sideplate,Fixer(DS-M6 series)	\$26.96
473549*			
286904	T286904	Sideplate,Developer(NDS-M6 series)	\$29.14
473550*			
286907	T286907	Sideplate, Developer (DS-M6 series)	\$33.31
473552*			
286928	T286928	Sideplate, Wash (NDS-M6 series)	\$44.16
473553*			
286930	T286930	Sideplate,Wash(DS-M6 series)	\$68.67
473554*			
286914	T286914	Sideplate, Fixer (NDS-M6 series)	\$27.30
473993*			
246879	T246879	Sideplate, Entrance	\$8.94
474100	T474100	Spacer	\$10.00
474114*			
990795	T990795	Isolation Relay	\$669.55
474115*			
990795	T990795	Isolation Relay	\$669.55
474118	T474118	Isolation Relay Harness	\$94.94
474119	T474119	Switch Harness	\$127.89
474124	T474124	Bracket, Idle Side	\$186.42
474125	T474125	Bracket, Drive Side	\$312.54
474157	T474157	Front Name Plate (M6A series)	\$265.67
474180	T474180	Arm	\$11.33
474185	T474185	Block	\$24.23
474198	T474198	Blower Wheel	\$127.45
474323	T474323	Screw (each)	\$5.00
474345	T474345	Pad (M7B Recirc. Pump)	\$10.00
474362	T474362	Impeller	\$165.44
474388	T474388	Knurled Pulley	\$50.21
474413	T474413	M6N Belt	\$71.48
474442	T474442	Belt, Blower (M6 series)	\$18.70
474469*			
5B6854	T5B6854	"T" Roller	\$130.98
474472	T474472	Rod Assembly Guide	\$301.58

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
474478	R474478	Dev/Fix X-Over Assy. (M6A series-Rubber)	\$580.49
474478*			
N.A.	S474478	Dev/Fix X-Over Assy. (Rubber/RubberSpring-Loaded)	\$684.49
474481	R474481	Detector Assembly (M7B)	\$908.75
474483	R474483	Developer Rack (M7B/M20)	\$2,881.09
N.A.	S474483	"Narrow Chain" Mammo Dev. Rack (M7B/M20)	\$4,685.00
474596	T474596	Bracket, Guide Shoe	\$10.00
474599	T474599	"N" Rewet Roller	\$218.17
474619	T474619	Roller-(Silicone-Wash Rack-M6A/M6B/460-480RA)	\$457.76
474621	T474621	Roller-(Silicone-Wash Rack-M6A/M6B/460-480/RA)	\$493.61
474645	T474645	Spring	\$10.00
474654	T474654	Impeller	\$91.92
474692	T474692	Pulley, Blower	\$126.76
474693	T474693	Pulley, Dryer Motor	\$276.32
474694	T474694	Blower Motor (Pulley) Assy.	\$657.60
474777	T474777	Stud-Snap	\$5.00
474836	T474836	Tubing, Anode Fix (M6B series)	\$24.16
474837	T474837	Tube-Anode Fix Htr.	\$21.02
474947	R474947	Reconditioned* M7A PCB 100 Board	\$623.83
474948*			
696164	T696164	Developer Heater (M7 series)-1000W	\$552.51
475409	T475409	Circuit Breaker (CB2/3)	\$87.19
475453	T475453	"T" Roller (Silicone)	\$219.02
475492	T475492	Cap ³ / ₄ Vinyl	\$5.00
475607*			
566313	T566313	Receptacle, TDR (M6 series)	\$10.07
475931	T475931	Overtemp Thermostat Box Assy.	\$98.66
476330	T476330	Fuse Holder	\$10.00
476350*			
5E1878	T5E1878	Switch-1	\$238.64
476378	T476378	Capacitor (ML series)	\$24.10
476462	T476462	Breaker, CB1 (M7B)	\$236.69
477000	T477000	Sideplate, Non-Drive	\$10.00
477001	T477001	Sideplate, Drive Side	\$10.00
477005	T477005	Sideplate, X-Over (M7,M35 series)	\$10.00
477011	T477011	Y Connector (1/2)	\$18.29
477012*			
914429	T914429	Replen. Pump-230V (M6AN/M6B series)	\$2,303.94
477017	T477017	Cover	\$49.64
477029	T477029	Stop	\$5.00
477062	T477062	Magnetic Detector Switch Kit	\$827.62
477092	T477092	Spacer, Split	\$5.00
477093	T477093	Spacer	\$10.00
477142	R477142	Reconditioned Heat Exchanger Guard	\$213.00
478548*			
990419	T990419	Nut, Tank Bulk, lrg.	\$10.00
478549*			
990400	T990400	Fitting	\$9.15
478560	T478560	Switch-PAP CLR, JOG	\$31.74
479788	T479788	Tube-Recirc. Pump	\$8.46
479815	T479815	Elbow 45° Angle	\$24.16

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
479820	T479820	Spring, Roller	\$12.22
479821	T479821	Graduated Cylinder	\$83.86
480059	T480059	Toggle Switch	\$32.14
480113	T480113	Bonnet-Drain Valve	\$33.77
480131	T480131	Motor Fan (M35 series)	\$31.82
480168	T480168	Chain, Main Drive (M35 series)	\$31.50
482222	T482222	Shaft, Drive Main (M35 series)	\$309.14
482271	T482271	Washer, Plain	\$5.00
482272	T482272	Spacer, Idler Gear	\$10.00
482273	T482273	Collar, Dryer Drive Shaft	\$63.63
482274	T482274	Spacer-Upper Sq. Entrance	\$10.00
482275	T482275	Gear, Dryer Drive	\$63.63
482277	T482277	Spacer-Up Sq. Exit	\$10.00
482283	T482283	Nut, Special (M35 series)	\$7.85
482339	T482339	Clip-Tank End (M7B)	\$45.49
482362*			
3H4475	T3H4475	Main Drive Motor-60hz (M6A series)	\$1,098.29
482375	T482375	Capacitor Guard Assy.	\$243.24
483007	T483007	Tape-Foam-Gray (per foot)	\$5.00
485291	T485291	Handle-Support	\$34.02
485391	T485391	Dryer Cover, Receiving End (M35 series)	\$128.08
485399	T485399	Heater (120V)	\$417.12
485406	T485406	Silicone Cross-Over Roller (M35 series)	\$307.94
486265	T486265	Screw (CB1)	\$5.00
486472	T486472	Screw	\$10.00
486480	T486480	Elbow	\$5.00
488163*			
651861	T651861	Replen. Pump-230V (M35A/M)	\$2,099.79
488169	T488169	Spacer (M35 series)	\$10.00
488190	T488190	Formed Tubing	\$15.49
489163	T489163	Shaft	\$51.03
489174	T489174	Washer, Rubber	\$25.25
489183	T489183	Insert-Bushing (M6A series)	\$3.68
489188	T489188	Stud Removal Tool	\$13.51
489693*			
1C7381	T1C7381	Lustran Utility Roller	\$36.98
489694*			
1C7387	T1C7387	Lustran "H" Roller	\$23.34
489695*			
1C7387	T1C7387	Lustran "H" Roller	\$23.34
489696*			
1C7390	T1C7390	Lustran Roller	\$24.23
489697*			
1C7389	T1C7389	Lustran "E" Roller	\$23.82
489698*			
1C7383	T1C7383	Roller (M7 series)	\$55.57
489699*			
1C7390	T1C7390	Lustran Roller	\$24.23
489700	T489700	Roller (M6 series)	\$33.90
489701*	T264425	Exhaust Hose (M35 series)	\$31.69

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
264425			
489702*			
1C7396	T1C7396	Roller	\$33.78
489703*			
1C7397	T1C7397	"W" Roller	\$31.82
489721	T489721	Shim, Dev. Tank	\$12.08
489788	T489788	Formed Tubing	\$9.56
489790*			
990417	T990417	Nut, Tank Fitting	\$10.72
489791*			
990418	T990418	Nut, Gland	\$10.00
489792*			
990400	T990400	Fitting	\$9.15
489793*			
990419	T990419	Nut, Tank Bulk, Fitting	\$10.00
489794*			
990415	T990415	Fitting-Recirc. Small	\$9.49
489795*			
990416	T990416	Fitting, Recirc., lrg.	\$13.10
489797	T489797	Protector-Feed Shelf (M20)	\$46.31
489997	T489997	Standoff, Detector Assembly	\$20.00
490092		Bearing, D/F X-Over (M35 series)	\$12.35
490096	T490096	Screw	\$5.00
490097	T490097	Screw	\$5.00
490108	T490108	Roller, Entrance Assy.	\$201.98
490109	T490109	Roller, Entrance Assy.	\$269.33
490160*			
9B8719	T9B8719	Tan Bearing	\$6.83
or	S779119		\$10.97
492100*			
1C7387	T1C7387	Lustran "H" Roller	\$23.34
492101*			
1C7381	T1C7381	Lustran Roller	\$36.98
492102*			
1C7390	T1C7390	Lustran Roller	\$24.23
492103*			
1C7389	T1C7389	Lustran "E" Roller	\$23.82
492104*			
1C7386	T1C7386	Lustran "D" Roller	\$64.01
492227*			
914429	T914429	Replenishment Pump (M6AN/AW; M6B series)	\$2,303.94
492260	T492260	Seal – Heat Exchanger	\$5.00
492948	T492948	Resistor	\$50.27
493505	T493505	Washer-Thrust	\$6.89
493674	T493674	Elbow (3/8x5/8)	\$10.00
495060	T495060	Cable Assy.-Magnetic Detector (M7B)	\$19.79
495071	T495071	Thermistor (10K)	\$318.15
495748	R495748	Fixer Rack (M7B/M20)	\$4,883.05
495749	R495749	Wash Rack (M7B)	\$2,900.21
499753	T499753	Retainer (3/16)	\$5.00
499754	T499754	Hose Connector	\$5.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
5B6001	T5B6001	Guide Shoe Assy. (270RA)	\$207.14
5B6058	T5B6058	Roller Assy., Stainless Steel "T"	\$570.65
5B6061	T5B6061	Roller, Stainless Steel, Dryer (180LP)	\$231.15
5B6064	T5B6064	Roller, Stainless Steel (270RA)	\$591.32
5B6076	T5B6076	"B" Roller (phenolic) (M6 series)	\$252.44
5B6077	T5B6077	"B" Roller (phenolic) (M35 series)	\$233.98
5B6099	T5B6099	Roller, Silicone (5000RA)	\$582.41
5B6151	T5B6151	Tubing, Formed	\$31.82
5B6163	R5B6163	Detector Entrance Assy. (M35A series)	\$2,164.76
5B6165	T5B6165	Extension Spring	\$11.06
5B6300	T5B6300	Tank Fitting	\$123.17
5B6301	T5B6301	Filter Tube	\$16.11
5B6316	T5B6316	Screw	\$5.00
5B6352	T5B6352	Recirculation Pump (5000RA)	\$1,192.77
5B6370	T5B6370	Roller	\$91.04
5B6371	T5B6371	Shaft, Roller	\$111.51
5B6395	T5B6395	Pin w/flat (3000RA)	\$10.00
5B6396	T5B6396	Pin-Shortened	\$10.00
5B6406	T5B6406	Roller, Stainless Steel "F" (3000RA)	\$591.32
5B6407	T5B6407	Roller, Rubber "G" (3000RA)	\$161.03
5B6409	T5B6409	Guide Shoe, Turnaround-Side	\$964.60
5B6411	T5B6411	Guide Shoe, Turnaround-Bottom	\$124.49
5B6414	T5B6414	Roller, Stainless Steel (3000RA)	\$606.31
5B6415	T5B6415	Buffed "D" Roller (3000RA)	\$168.34
5B6527	T5B6527	Rocker Bearing	\$12.29
5B6569	T5B6569	Elbow	\$10.00
5B6575	R5B6575	Wash Rack (3000RA)	\$5,334.66
5B6576	R5B6576	Dev/Fix X-Over Assy. (3000RA)	\$1,321.61
5B6587	T5B6587	Shaft, Idle (3000RA)	\$43.85
5B6588	T5B6588	Shaft, Driven (3000RA)	\$54.56
5B6589	T5B6589	Spacer	\$32.15
5B6590	T5B6590	Roller	\$112.33
5B6591	T5B6591	Gear, Wash (3000RA)	\$14.47
5B6630	T5B6630	Dole Valve (.50 gpm)	\$128.27
5B6682	T5B6682	Solenoid, Water Inlet	\$190.64
5B6775	T5B6775	Spacer	\$22.59
5B6776	T5B6776	Spacer	\$84.92
5B6804	T5B6804	Tube	\$27.57
5B6808	T5B6808	Shoulder Screw (3000RA)	\$13.46
5B6817	T5B6817	Shaft, Roller Idler (3000RA)	\$143.39
5B6850*			
166331	T166331	Nut	\$5.00
5B6854	T5B6854	Roller, Stainless Steel "T" (with gudgeons)	\$130.98
5B6890	R5B6890	Wash Rack Assy. (270RA)	\$4,547.89
5B6913	T5B6913	Formed Tube	\$162.29
5B6917	T5B6917	Formed Tube	\$72.77
5B6929	T5B6929	Tube Assy. (5000RA)	\$470.86
5E1878	T5E1878	Switch	\$238.64
5E2132	T5E2132	Spacer, Master Roller (outer) (Min-R)	\$9.83
5E2133	T5E2133	Spacer, Master Roller (inner) (Min-R)	\$95.45
500340*			

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
616464	T616464	Circuit Breaker-2 amp	\$39.75
502464	T502464	Pin	\$5.00
508498*			
479820	T479820	Spring, Roller	\$12.22
512893	R512893	Wash Rack Assy. (M6AN)	\$1,777.83
516425	T516425	Plate, Floor (M7B)	\$13.65
520388	T520388	O-Ring	\$17.20
520389*			
520391	T520391	Nut	\$31.82
520394	T520394	Washer	\$10.00
521417	T521417	Bulb (for ID camera)	\$263.40
522368*			
7C6650	T7C6650	Recirc. Pump Guard	\$74.28
527068	T527068	Tube	\$10.00
529988	T529988	Boot-Plug	\$38.62
531263	T531263	Orifice, T-Tray	\$38.18
531270	T531270	Tube	\$16.18
532007	T532007	Valve Body	\$39.69
532014*			
911809	T911809	Sideplate, Fixer (DS-M7 series)	\$25.32
532015*			
911845	T911845	Sideplate, Fixer (NDS-M7 series)	\$25.25
532016*			
286933	T286933	Sideplate, Developer (DS-M7 series)	\$43.41
532017*			
286934	T286934	Sideplate, Developer (NDS-M7 series)	\$94.12
532018*			
286924	T286924	Sideplate, Wash (DS-M7 series)	\$26.89
532019*			
286926	T286926	Sideplate, Wash (NDS-M7 series)	\$27.71
532033	T532033	Adapter	\$7.00
532037	T532037	Roller	\$94.50
532054	T532054	Bearing (M35, Dryer)	\$7.78
532062	T532062	Chain-w/14mm-wide links (M7 Dev.)	\$60.48
532070*			
SP532075	T532075	O-Ring	\$10.00
532086	T532086	Detector Roller, Lower	\$390.03
532087	T532087	Detector Roller, Upper	\$350.78
532092	T532092	Clip, Hold Down	\$22.59
532095	T532095	Heat Exchanger (M6 series)	\$103.45
532105	T532105	Grommet (Laser)	\$1.18
532124	T532124	Gear-Snap (22-tooth)	\$8.87
532155	T532155	Roller Assy.	\$660.89
532159	T532159	Switch Assembly	\$46.81
532164	T532164	Tie Rod	\$292.70
532171	T532171	Timing Belt (M7 series)	\$10.71
532176	T532176	Bearing, Dryer (M7 series)	\$5.00
532180	T532180	Timing Belt (M7 series)	\$23.14
532181*			
489163	T489163	Shaft	\$51.03
532204	T532204	Bearing	\$155.30

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
532219	T532219	Shim	\$205.57
532229	T532229	Gear Sprocket (M7B)	\$95.13
532247	T532247	Sprocket, 19-tooth-50 Hz (M7B)	\$448.50
532248	T532248	Chain, Main Drive (M7 series)	\$95.45
532255	T532255	Spacer, Defl. Refl.	\$31.82
532289	T532289	Shaft-Dry Drive (M7B)	\$269.64
532290	T532290	Roller Assy.-Dryer Ent. (M7B dryer)	\$245.89
532301	T532301	Spacer, Drive G (M7 series)	\$18.29
532302	T532302	Gear, Dry Drive (M7 series)	\$61.55
532337*			
594762	R594762	Reconditioned Blower (M7B)	\$1,458.19
532361	T532361	Roller Assy. (Dryer)	\$115.23
532363	T532363	Sprocket, Main Drive (M7B)	\$60.29
532378	T532378	Seal (O-Ring)	\$4.64
532379	T532379	Washer (each)	\$5.00
532739	T532739	Timer, Standby (A-1, 30 minute)	\$617.21
534714	T534714	O-Ring	\$5.00
534738	T534738	Tank Clip (M7B)	\$318.15
534739*			
482339	T482339	Clip – Tank End (M7B)	\$45.49
534740*			
482339	T482339	Clip – Tank End (M7B)	\$45.49
534742	T534742	Tank End Seal ("T")	\$10.00
534753	T534753	Shim – Weir	\$11.40
534789*			
650710		Bearing	\$10.00
534790*			
477092	T477092	Spacer, Split	\$5.00
534791	T534791	Spacer, Rack D.S.	\$11.06
535004	T535004	Lockwasher	\$5.00
535639	T535639	O-Ring	\$9.21
536834	T536834	K-6 Relay Kit (M7A)	\$79.10
536837	T536837	Dryer Knob (M7B)	\$98.66
536842*			
477093	T477093	Spacer	\$10.00
536855	T536855	M7B Switch Interlock Assy.	\$114.35
536884	T536884	Rubber "B" Roller	\$313.36
539110	T539110	Clip	\$33.58
539134	T539134	Tube, Formed (M7 series)	\$42.90
539138	T539138	O-Ring (Weir side)	\$5.00
539139	T539139	Tubing (12.25")	\$14.81
539141	T539141	Tubing	\$16.52
539151	T539151	Orifice (M35 series)	\$14.13
539152	T539152	Orifice (M7B)	\$99.41
539197	T539197	O-Ring, Def. Roll.	\$5.00
539198	T539198	Gudgeon, Plain	\$9.28
541891	T541891	Plug, thermowell	\$5.00
542092	T542092	Cover Stop Strap (M7B)	\$27.50
543884	T543884	Pin	\$24.16
543993	T543993	Bearing (M35 Dryer)	\$5.00
543995	T543995	Roller (M20 Dryer)	\$258.17

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
544038*			
532033	T532033	Adapter	\$7.00
544081*			
SP164385	T164385	Pin	\$10.00
544127	T544127	Drain Valve (<i>fix. & dev.</i>)	\$94.25
544154	T544154	Lamp	\$29.62
544260	T544260	Elbow (<i>M35 series</i>)	\$10.00
544301	T544301	Pin Housing	\$10.00
544302	T544302	Plug	\$10.00
548952	T548952	Spacer, Deflector Roller	\$20.68
548995	T548995	Screw	\$5.00
550774	T550774	Gear. Drive Rack (<i>Black</i>)	\$7.03
550781	T550781	Drive Gear (<i>White</i>)	\$10.03
551400	T551400	Elbow	\$10.00
554780		Seal, Tank End ("L")	\$10.00
556104	T556104	Adapter	\$5.00
556289	T556289	Nut	\$5.00
557422	T557422	PCB 100-M7B, M35/M35A	\$2,142.53
558655	T558655	Connector (<i>M6B series</i>)	\$27.37
558690	T558690	Connector, P29, P26, P21, P17	\$7.94
559569	T559569	Cable Clamp	\$12.55
559787	T559787	Shaft, "G" Roller	\$51.91
562244	T562244	Buzzer DS5	\$85.74
562289	T562289	Filter Line, FL1, FL2	\$246.90
562334	T562334	Relay, K2-K8 (<i>M7B series</i>)	\$127.89
562365	T562365	Switch, Water Fl.	\$509.67
562366	T562366	Fix/Dev Thermistor Assy. (<i>M6B series</i>)	\$532.54
562367	T562367	Overtemp Thermistor (<i>M6B series</i>)	\$322.18
562389	T562389	Guide Shoe, D/F X-Over	\$485.79
562418	T562418	Suppression-Network	\$15.70
562477*			
594394	T594394	Replen. Pump, M8 (240V)	\$1,692.42
562485	T562485	O-Ring, Filter Seal (<i>M6B series</i>)	\$10.37
562500	T562500	Spring-Extension	\$114.53
562602*			
990418	T990418	Nut, Gland	\$10.00
562613	T562613	Spacer	\$20.68
562614	T562614	Shaft (<i>270RA series</i>)	\$94.25
562619	T562619	Chain-26 link	\$15.42
562641	T562641	Spacer	\$31.82
562693	T562693	Relay (K1)	\$66.97
562790	T562790	Developer Filter (<i>M6 and M7 series</i>)	\$10.00
562810	T562810	Solenoid Valve (<i>M8 series</i>)	\$149.56
562864	T562864	Shaft, Squeegee (<i>270RA series</i>)	\$548.54
562891*			
990418	T990418	Nut, Gland	\$10.00
562892*			
990419	T990419	Nut-Tank Bulk Fitting large	\$10.00
562913	T562913	Rubber Cover Roller	\$92.61
562956	T562956	ID Camera Board (<i>Serial Number 7000 and above</i>)	\$721.93
563058	T563058	Clamp	\$5.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
563689	T563689	Dryer Roller, Rubber (<i>M7B series</i>)	\$247.46
566313	T566313	Receptacle, TDR (<i>M6 series</i>)	\$10.07
570450	T570450	Back-up Plate	\$10.00
571241	T571241	Bearing, Dryer (<i>M7B/M35 series</i>)	\$10.00
573971	T573971	O-Ring	\$5.00
575801	T575801	Bearing Rocker	\$12.83
575812	T575812	Detector Switch Assy.	\$71.25
575820	T575820	Screw	\$5.00
575821	T575821	Screw	\$5.00
575822	T575822	Washer	\$5.00
575834	T575834	Standoff, S.S.	\$6.83
575903	T575903	Timer, 3-minute	\$258.24
579383	T579383	Fitting – Elbow	\$10.00
579424	T579424	Circuit Breaker (CB3)	\$111.01
579427	T579427	Circuit Breaker-3 amp (CB2)	\$94.37
581200	T581200	Drain Valve	\$92.67
581241	T581241	Gasket	\$10.00
581243*			
492260	T492260	Seal Heat Exchanger	\$5.00
583127	T583127	Valve, Chemical Drain	\$106.79
583177	T583177	Varistor-RV3 (<i>TDR</i>)	\$13.09
583178	T583178	Varistor-RV1 (<i>Trans. Term.</i>)	\$25.94
583514	T583514	Seal, Tank Separator (<i>M7 series</i>)	\$10.00
585721	T585721	Knob (<i>M6B series</i>)	\$29.36
585758	T585758	Terminal	\$5.00
587329	T587329	Dryer Roller (<i>M7 series</i>)	\$99.67
591139	T591139	Filter, Air	\$272.92
592372	T592372	Formed Tubing (<i>M7 series</i>)	\$21.23
592373	T592373	Formed Tubing (<i>M7 series</i>)	\$15.36
592374	T592374	Tube (<i>M7 series</i>)	\$41.71
N/A	T593739	Timer It's Back!	\$181.94
593740*			
1C7772	T1C7772	Relay	\$327.79
594229	T594229	Shaft, Drive	\$36.54
594230	T594230	Shaft, Driven	\$53.30
594234	T594234	Roller Assy., Deflector	\$60.80
594236	T594236	Shaft, Fixer (<i>270RA series</i>)	\$35.28
594237	T594237	Roller, Crossover (<i>270RA</i>)	\$97.34
594242	T594242	Shaft, Master Roller (<i>270RA</i>)	\$83.22
594254	T594254	Gear (20 tooth)	\$47.75
594264	T594264	Chain – 70 Link Developer (<i>M8</i>)	\$68.54
594265	T594265	Chain – 55 Link Fixer (<i>M8</i>)	\$38.45
594266	T594266	Chain – 42 Link (<i>M8</i>)	\$81.92
594271	T594271	Rewet Roller (<i>M8</i>)	\$114.60
594288*			
1C7551	T1C7551	Lustran "D" Roller	\$247.91
594293	T594293	Shaft-Rack Driver (<i>270RA</i>)	\$150.95
594294	T594294	Spacer	\$17.20
594388	T594388	Valve-Drain	\$155.99
594391	T594391	Recirculation Pump (<i>M8</i>)	\$3,256.86

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
594394	T594394	Replenishment Pump (M8)	\$1,692.42
594450	T594450	Potentiometer, Dryer	\$32.07
594453*			
7E6924	T7E6924	Temp. Display Module (M6B series)	\$1,370.70
594470	T594470	Switch	\$21.71
594479*			
748955	T748955	Worm Gear, L.H.	\$10.00
594480*			
748956	T748956	Worm Gear, R.H.	\$10.00
594534	T594534	Washer Coupler-Tank	\$10.00
594535	T594535	Gasket Coupler-Tank	\$5.00
594752	T594752	Terminal Strip (TB2)	\$88.45
594759*			
1C0033	T1C0033	Buzzer Assembly (M7 series)	\$35.41
594760	T594760	Replen. Pump-115V (M7 series)	\$2,041.46
594762	R594762	Reconditioned Blower (M7 series)	\$1,458.19
594764*			
594765	T594765	Developer Heater (M7 series)	\$579.75
594777*			
246868	T246868	Relay Kit (M7B, M35, M35A series)	\$415.49
594793*			
594796*			
594760	T594760	Replen. Pump-115V (M7 series)	\$2,041.46
594809*			
8B7075	R8B7075	Reconditioned Dryer Rack Assy. (M7B)	\$5,628.20
595426*			
192019	T192019	Bearing, Large	\$53.80
595428	T595428	Formed Tube (M7 series)	\$26.41
595860	T595860	Knob, Lock Control (M6B series)	\$52.48
596604	T596604	Tubing (Heat Exchanger)	\$7.99
596611	T596611	Duct, Flexible-Heat Exchanger (M6AW)	\$293.71
596646	T596646	Fan (M6AW)	\$313.36
597064*			
286910	T286910	Sideplate, Wash NDS-M6AW)	\$25.46
597065*			
286912	T286912	Sideplate, Wash (DS-M6AW)	\$24.91
597075	T597075	Formed Tube (M6 series)	\$9.56
597078	T597078	Tube (36 in.)	\$7.78
597096	T597096	Spacer, Split	\$27.57
597360	T597360	Timer (M6AW standby)	\$373.65
598815	T598815	Housing-Pin	\$10.00
N.A.	T598815P	Pin-Male	\$10.00
598816	T598816	AMP Connector-Male	\$5.00
N.A.	T598816P	Pin-Female	\$2.21
598828	T598828	Quencharc	\$9.21
598833	T598833	Grommet (Blower)	\$5.00
598834	T598834	Latch, Cabinet (M7 series)	\$7.03
598835	T598835	Collar (M7B)	\$15.56
598847	T598847	Dev. Thermostat-Probe (M7B)	\$548.48
603510*			
1F5313	T1F5313	Relay and Varistor	\$228.75

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
607696	T607696	Formed Tubing-Flow(M7 series)	\$23.07
607833	T607833	Gudgeon Sprocket	\$6.89
609952	T609952	Flow Control Valve (1.5 gpm)	\$142.19
611484	T611484	Thermistor Assy. M7A)	\$172.63
611486	T611486	ID Camera Board (Serial Number 6999 and lower)	\$4,633.34
611497	T611497	Varistor-Suppressor	\$10.00
611498	T611498	Varistor	\$8.74
613114	T613114	Gasket, Replenisher Filter	\$5.00
613122	T613122	Cover, Evaporation	\$672.21
613342	T613342	Stem & Washer Assy.	\$46.12
616122	T616122	Receptacle, AC Input	\$7.37
613276*			
990415	T990415	Fitting-Recirc. Small	\$9.49
613277*			
990416	T990416	Fitting-Recirc. Large	\$13.10
613278*			
990417	T990417	Nut-Tank Fitting	\$10.72
616289	T616289	Power Cord (I.D. Camera)	\$63.13
616398	T616398	Switch, S6 (M6B/M7B series)	\$43.28
616463	T616463	Circuit Breaker-15 amp	\$39.75
616464	T616464	Circuit Breaker-2 amp	\$39.75
618502	T618502	Diode	\$31.53
619365	T619365	Washer	\$5.00
619512	T619512	Washer	\$10.00
621230	T621230	Flow Indicator Assembly (sight tube)	\$41.27
621233	T621233	Water Solenoid (M7B)	\$357.46
621235	T621235	Weep-Line Orifice (M7 wash)	\$218.74
621472	T621472	100 Board-Developer Thermostat	\$389.72
621480*			
651125	T651125	Belt, Dryer Transport (M6 series)	\$80.01
622017	T622017	Recirc. Pump-230V (M6B series)	\$2,015.82
622018*			
914429	T914429	Replen. Pump-230V (M6AN/AW)	\$2,303.94
622019*			
748683	T748683	Blower Motor Assembly	\$838.07
622020	T622020	Buzzer (M6B series)	\$26.89
622022	T622022	DC Harness Assy.-Elect. Box (M6B series)	\$339.82
622024	T622024	Cable Assy.-Aux./Run (M6B series)	\$68.73
622025	T622025	Cable Assy.-Female to T622024 (M6B series)	\$67.42
622030	T622030	Heater Assy.-Developer	\$999.88
622036	T622036	Dryer Heater Wires	\$379.39
622040	T622040	Thermistor Assy. - Dryer (M6B series)	\$147.23
622046	T622046	ED Cable Assy.-Dryer Therm.	\$109.62
622049	T622049	ED-Cable Det. Sw. W/P (M6B series)	\$23.48
650604	T650604	Compressor Repair Kit (Laser)	\$309.83
650628	T650628	Tube, Replenishment Pump	\$187.05
650663	T650663	Gear, Drive	\$9.49
650672*			
748683	T748683	Blower Motor Assembly	\$838.07
650705	T650705	Bearing	\$11.06

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
650707	T650707	Bearing	\$10.00
650710		Bearing	\$10.00
650711	T650711	Bearing	\$11.12
650733*			
532219	T532219	Shim	\$205.57
650852	R650852	<i>Reconditioned</i> Developer Thermostat	\$2,000.17
650852	R650852-w/o	<i>Reconditioned</i> Developer Thermostat (<i>without thermistor</i>)	\$1,662.07
N.A.	R650880	<i>Reconditioned</i> Dev. Controller Assy. (M20)	\$1,267.92
650938	T650938	Gasket (2 ea.-10 ft. rolls)	\$27.05
651107	T651107	E-Ring	\$5.87
651125	T651125	Belt, Dryer Transport (M6 series)	\$80.01
651853	T651853	Buzzer (M35 series)	\$112.46
651860	T651860	Run/Standby Switch	\$172.43
651861	T651861	Replenishment Pump (230 volt, M35A/M20)	\$2,099.79
651862	T651862	Heater -220V (M35 series/2000)	\$419.52
651868	T651868	O.T. Thermostat (M35A series)	\$602.28
651869	T651869	Thermistor, 10K	\$150.70
651871	T651871	ED Wire Dryer Control	\$10.24
651872*			
779524	T779524	Potentiometer (M35 series)	\$583.33
651873	T651873	Quencharc-(pkg. 2)	\$201.66
651881	T651881	Door Chain (M35 series)	\$10.00
651893	T651893	Main Wire Harness Assy. (M35A series)	\$429.60
652486	T652486	Counter, (M6B)	\$187.32
655290*			
7B8547*			
1F5313	T1F5313	Relay, Solid State	\$228.75
655658	T655658	M6B Main Drive Motor-(S/N 9,999 & below)	\$1,250.46
655688	T655688	Tubing, Formed (M6B series) (18.5")	\$10.72
655700	T655700	M6B Valve-Dual Solenoid Inlet	\$98.28
659546	T659546	Supply-Power (M6B series)	\$714.43
N.A.	T66361-2	Pin Contact for Amp Plugs	\$1.09
696088	T696088	Circuit Breaker (2000A series)	\$123.35
696092	R696092	Exit Rack (270RA stand-alone)	\$4,643.86
696164	T696164	Heater, Solution (1000W) (M6B series)	\$552.51
696165*			
3H4056	T3H4056	Developer Heater (270RA)	\$720.20
696174*			
968802	T968802	800 Circuit Board & Cable Assy. (M35 series)	\$622.89
699901	T699901	Tube-Dev. Drain	\$8.33
7B7229	T7B8229	Captive Screw	\$7.85
7B8547*			
1F5313	T1F5313	Relay, Solid State	\$228.75
7C1222	T7C1222	Mount, Isolation (Laser)	\$26.48
7C6650	T7C6650	Recirc. Pump Guard	\$74.28
7C6660	T7C6660	Tie Rod (M35 series)	\$25.59
7C6661	T7C6661	Tie Rod ((M6 series)	\$32.56
7C6662	T7C6662	Tie Rod	\$24.50
7C6980	T7C6980	Replenishment Timer (M35 series)	\$439.05
7C7860	T7C7860	Guideshoe, Smooth (M35A series)	\$333.90
7C7862	R7C7862	Developer Rack (M35M/M35A-M series)	\$3,009.35

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
7C8226	T7C8226	Stand-Off (for 8B6902 thermostat)	\$5.00
7C8328	R7C8328	Developer Rack (3000RA)	\$8,089.70
7C8476	T7C8476	Roller, S.S. Dev. Rack (3000RA)	\$607.95
7C8477	T7C8477	Roller, Dev. Rack (3000RA)	\$240.16
7C8478	T7C8478	Water Solenoid (M6 series)	\$3,134.32
7E6924	T7E6924	Temp. Display Module (M6B series)	\$1,370.70
7E6936	S7E6936	Dryer O.T. (for TI-BA TBR09 ONLY!**)	\$35.03
748021	T748021	Shaft, Gear Dryer (270RA)	\$134.63
748022	T748022	Shaft, Idle Dryer (270RA stand-alone/integrated)	\$112.33
748126*			
1C7558	T1C7558	Lustran "A" Roller (M35 series)	\$6105
748138	T748138	Stiff-Rib (M6 series)	\$201.35
748147	T748147	Roller, X-Over (phenolic) (M35/M series)	\$121.91
748164	R748164	Turnaround Assy. (M35 series)	\$689.88
748165	R748165	Fixer Rack Assembly (M35A series)	\$3,755.71
748245	T748245	Clamp	\$32.01
748248	T748248	Clamp (M6B series)	\$68.17
748249	T748249	Elbow	\$8.94
748426	T748426	Access Panel (M6B series)	\$816.59
748427	T748427	Nipple	\$10.00
748430	T748430	Tubing, Formed (M6B series)	\$9.56
748437	T748437	Dryer Pan	\$283.25
748448	T748448	Tubing, Formed	\$8.53
748511	T748511	Heater Element (270RA)	\$880.75
748585	T748585	Connector Tee	\$11.67
748589	T748589	Hose, Formed (M6B series)	\$25.73
748683	T748683	Blower Motor Assembly with cord	\$838.07
N.A.	S748683	Blower Motor Assembly (50 hz-European Version)	\$838.07
748954	T748954	Drive Shaft (M7B)	\$329.43
748955	T748955	Worm Gear, L.H.	\$10.00
748956	T748956	Worm Gear, R.H.	\$10.00
748957	T748957	Washer-Thrust (M6/M7 series)	\$5.00
752731	T752731	Drive Chain (Min-R series)	\$31.53
760476	T760476	Formed Tubing (1/2 I.D.)-per foot	\$5.00
760477	T760477	Formed Tubing (3/4 I.D.)-per foot	\$9.76
762095*			
760476	T760476	Tubing (1/2") (per foot)	\$5.00
N.A.	T763106	Paint (9-oz. Spray), Sand Beige	\$57.14
N.A.	T763107	Paint (9-oz. Spray), Image Light	\$57.14
N.A.	T763108	Paint (9 oz. spray), Greige	\$57.14
764106	T764106	Terminal (pkg. of 10)	\$10.00
764111	T764111	Terminal(pkg. of 10)	\$11.12
764390	T764390	M6AW Pin (pkg. of 10)	\$10.00
779119*			
9B8719	S779119	"Special" (TGM-gray) Bearing	\$10.97
779119*			
9B8719	T9B8719	Tan Bearing	\$6.83
779120	T779120	Drive Shaft & Worm Assy (M7 series)	\$601.78
N.A.	T779120C	Complete Main Drive Assy. (M7 series)	\$660.00
779155	R779155	Fix/Wash X-Over Assy. (M6A series)	\$1,004.73
779188	R779188	Fix/Wash X-Over Assy. (460RA,480RA,M6I)	\$1,177.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
779216	T779216	Rack Bar (M35 series)	\$90.59
779241	T779241	Valve, Water Inlet	\$88.70
779243	T779243	Valve, Water Dump (270RA/3000)	\$210.61
779332*			
IC4593	T1C4593	Replenishment Pump (180,270,460,480,3000,5000RA)	\$1,808.15
779334	T779334	Chain	\$42.65
779344	T779344	Water Solenoid Valve	\$990.99
779349	T779349	Splash Guard, Wash (270RA/180LP/3000RA)	\$318.15
779360	T779360	Tubing, Formed	\$24.91
779362	T779362	Tubing, Formed	\$40.64
779363	T779363	Tubing, Formed	\$20.68
779364	T779364	Tubing, Formed	\$12.29
779365	T779365	Tubing, Formed	\$15.49
779366	T779366	Tubing, Formed	\$17.20
779368	T779368	Tubing, Formed	\$17.20
779369	T779369	Tube	\$11.67
779370	T779370	Tube	\$10.10
779371	T779371	Tube	\$16.18
779382	T779382	Formed Tube	\$8.67
779383	T779383	Formed Tube	\$13.79
779384	T779384	Formed Tube	\$20.68
779385	T779385	Tubing	\$10.37
779387	T779387	Tubing, Formed	\$20.68
779388	T779388	Tubing, Formed	\$37.67
779460	T779460	Tube	\$127.26
779461*			
943667	T943667	Tubing, Formed	\$20.68
779462	T779462	Tube	\$12.17
779501	T779501	Motor Capacitor (M35A/M)	\$186.42
779504	T779504	AC Motor Controller (M35A/M)	\$1,182.14
779506	T779506	Relay, K-1 (M35A/M)	\$65.84
779524	T779524	Potentiometer Assy. (M35A/M)	\$583.33
779562	T779562	Thrust Spacer (M7 drive)	\$30.99
779727	T779727	Pulley, Drive-Takeup	\$44.99
779728	T779728	Sleeve, Pulley	\$30.99
779729	T779729	Sprocket, Drive, 30-tooth	\$86.94
779749	T779749	Sleeve, Idler	\$32.47
788553	T788553	Stripping Foam (per foot)	\$32.83
8081176	T8081176	M35 Stand (with stationary shelf)	\$610.00
8261703	R8261703	Refurbished I.D. Camera, Model 2	\$2,992.76
8B6568	T8B6568	Squeegee "J" Roller	\$161.60
8B6577	R8B6577	Exit Rack (180LPS)	\$215.71
8B6583*			
8B6581	T8B6581	"D" Roller (M35)	\$239.59
8B6612	R8B6612	Wash Rack Mod Kits 29,19,7 (270RA)	\$2,189.19
8B6679	R8B6679	Fixer Rack Assy. (270RA)	\$6,695.47
8B6682	R8B6682	Detector X-Over Assy. (3000RA)	\$801.11
8B6844	T8B6844	Shaft	\$80.89
8B6953	T8B6953	Spacer (3000RA)	\$24.02
8B6955	T8B6955	Spacer	\$10.92
8B6965	R8B6965	Dev/Fix Cross-Over Assy. (5000RA)	\$2,931.04

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
8B6966*			
3H8726	R3H8726	F/W Cross-Over Assy. (5000RA)	\$3,566.93
8B6969	R8B6969	Wash Rack (5000RA)	\$4,592.28
8B6970	R8B6970	Squeegee Assy. (5000RA)	\$4,535.92
8B7050	T8B7050	Brkt.-Adj. (3000RA)	\$10.51
8B7060*			
1F5313	T1F5313	Solid State Relay (20 amp)	\$228.75
8B7064*			
1F5313	T1F5313	Solid State Relay (12 amp)	\$228.75
8B7075	R8B7075	Reconditioned Dryer Rack (M7 series)	\$5,628.20
8B7110	T8B7110	Squeegee Drive Roller	\$492.85
8B7111	T8B7111	Squeegee Roller Driven	\$488.63
8B7293	R8B7293	Developer Rack Assy. (5000RA)	\$4,883.30
8B7294	R8B7294	Fixer Rack Assy. (5000RA)	\$5,018.29
8B7405	T8B7405	Bearing (5000RA)	\$10.72
834692	T834692	Knob (M35 series)	\$47.19
SP851345	T851345	Fuse (250V)	\$5.00
852520	T852520	Screw	\$5.00
852522	T852522	Screw	\$5.00
852523	T852523	Screw	\$5.00
852525	T852525	Screw	\$5.00
852555	T852555	Set Screw	\$5.00
852558	T852558	Set Screw	\$5.00
852559	T852559	Set Screw	\$5.00
852611	T852611	Washer	\$5.00
SP852613	T852613	Washer	\$5.00
SP852615	T852615	Washer	\$5.00
852616	T852616	Washer	\$5.00
852617	T852617	Lock Washer	\$5.00
SP852621	T852621	Nut	\$10.00
SP852622	T852622	Nut	\$5.00
852623	T852623	Nut	\$5.00
SP852656	T852656	Lock Washer	\$5.00
852670	T852670	Hex Screw	\$5.00
SP852682	T852682	Screw	\$8.80
852772	T852772	Screw	\$5.00
852774	T852774	Screw	\$5.00
852775	T852775	Screw	\$5.00
852776	T852776	Screw	\$1.03
852779	T852779	Screw	\$5.00
852799	T852799	Washer	\$0.53
SP852801	T852801	Nut	\$5.00
852802	T852802	Screw	\$0.86
K852807	T852807	Screw	\$5.00
SP852808	T852808	Screw	\$5.00
852810	T852810	Screw	\$0.71
852811	T852811	Screw	\$5.00
852813	T852813	Screw	\$0.52
852821	T852821	Screw	\$5.00
852831	T852831	Washer	\$1.25
SP852835	T852835	Screw	\$5.00

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
852845	T852845	Set Screw	\$5.00
SP852847	T852847	Set Screw	\$5.00
852851	T852851	Screw	\$5.00
852868	T852868	Screw	\$5.00
852871	T852871	Screw	\$5.00
SP852885	T852885	Screw	\$5.00
SP852887	T852887	Screw	\$5.00
852888	T852888	Screw	\$0.88
852889	T852889	Set Screw	\$5.00
852895	T852895	Screw	\$5.00
852941	T852941	Screw	\$5.00
852946	T852946	Screw	\$0.32
SP852959	T852959	Screw	\$5.00
852960	T852960	Screw	\$0.36
852961	T852961	Lockwasher	\$5.00
SP852964	T852964	Screw	\$5.00
852971	T852971	Screw	\$5.00
852972	T852972	Screw	\$5.00
852973*			
SP165678	T165678	Standby Screw	\$5.00
852974	T852974	Screw	\$5.00
853000	T853000	Pin	\$5.00
SP853006	T853006	Screw	\$5.00
853032	T853032	Pin	\$5.00
SP853041	T853041	Screw	\$5.00
853051	T853051	Ring, Retaining	\$5.00
SP853069	T853069	Pin	\$30.78
853080	T853080	Hex Nut (M7B)	\$0.88
853087	T853087	Lock Washer	\$0.88
SP853096	T853096	Lock Washer	\$5.00
853098	T853098	Screw	\$0.69
SP853100	T853100	Pin	\$5.00
853141	T853141	Screw	\$0.88
853149	T853149	Set Screw	\$1.17
SP853165	T853165	Screw	\$5.00
853186	T853186	Hex Bolt	\$5.00
892099	T892099	Flow Control Valve	\$158.41
897205	T897205	Bearing-Thrust	\$9.11
9B8086	R9B8086	Dryer Rack (180LP)	\$3,674.52
9B8087	T9B8087	Squeegee Kit/Drive Roller	\$334.22
9B8088	T9B8088	Squeegee Kit/Driven Roller (M6 series)	\$467.71
9B8276	R9B8276	Developer Rack (270RA-integrated)	\$7,389.90
9B8293	T9B8293	Guideshoe, Smooth (M6 series)	\$489.57
9B8296	T9B8296	Adj. Fork	\$54.31
9B8359	T9B8359	Blower Motor (M6B series)	\$990.53
N.A.	S9B8359	Blower Motor (50 hz-European Version)	\$990.53
9B8436	R9B8436	Dryer Rack (270RA integrated)	\$7,588.94
9B8445	R9B8445	Fix/Wash X-Over Assy. (270RA)	\$2,347.42
9B8450	R9B8450	Fix/Wash X-Over Assy. (180LP)	\$766.31
9B8625	T9B8625	Spring	\$44.54
9B8626	T9B8626	Spring	\$80.20

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
9B8637	T9B8637	Roller, Stainless	\$577.79
9B8652*			
9B8436	R9B8436	Dryer Rack (270RA integrated)	\$7,588.94
9B8715	T9B8715	"D" Roller (3000RA)	\$574.25
9B8719	T9B8719	Tan Bearing	\$6.83
9B8720	T9B8720	Bearing	\$26.69
9B8828*			
3H8726	R3H8726	Fix/Wash X-Over (5000RA)	\$3,566.93
9B8832	R9B8832	Fix/Wash X-Over Assy. (3000RA)	\$2,698.48
9B8951	T9B8951	Chain (180, 270RA, 3000RA)	\$43.34
9B8960	T9B8960	Locking Plate	\$7.03
9B8964	T9B8964	Sprocket, Snap-On	\$15.90
9B8965	T9B8965	Sprocket, Drive	\$16.341
9B8967	T9B8967	Spacer-Main Drive	\$5.00
9B8968	T9B8968	Gudgeon, Chain Idler	\$31.26
911150	T911150	Drive Chain (480RA/5000RA)	\$34.02
911207	T911207	Nameplate (M7 series)	\$59.79
911218	T911218	Nameplate (M6B series)	\$92.23
911343	T911343	Mount-Isolator (Laser)	\$36.78
911405	T911405	Connector	\$17.20
911424	T911424	Filter Line	\$231.08
911483	T911483	Tube	\$6.89
911653	T911653	Shaft, Exit (3000RA)	\$97.65
911660	T911660	Roller Assy.-Phenolic (270RA)	\$89.78
911661*			
5B6058	T5B6058	Roller Assy. (stainless)	\$570.65
911725	T911725	Idler Bushing	\$12.83
911727	T911727	Drive Chain	\$50.09
911803*			
1C7554	T1C7554	"D" Roller (M35 series)	\$68.92
911804	T911804	Dev. Idler Roller (M35A/M)	\$184.53
911805	T911805	Roller Assy.-Cluster (M35A/M)	\$184.09
911806*			
1C7640	T1C7640	"B" Roller Rubber (M35 series)	\$319.91
911807	T911807	Roller (M35A/M)	\$298.05
911809	T911809	Sideplate, Fixer (DS-M7 series)	\$25.32
911845	T911845	Sideplate, Fixer (NDS-M7 series)	\$25.25
914160	T914160	Cotter Pin	\$5.00
914172	T914172	O-Ring	\$5.00
914187	T914187	Recirculation Pump (K-100)	\$1,095.11
914265	R914265	Dev. Turnaround Assy.	\$1,023.68
914267	T914267	Mammo Kit (M7,M7B,M35)	\$250.99
914268	T914268	Mammo Kit (M35/M35A/M20)	\$196.40
914310	T914310	Blower Bearing Guard (M6 series)	\$95.45
914429	T914429	Replen. Pump-230V (M6 series)	\$2,303.94
914441	T914441	Circuit Breaker- 20 amp (M7B)	\$195.05
914443	T914443	Circuit Breaker-30 amp (M6 series)	\$205.51
914444	T914444	Circuit Breaker-15 amp	\$185.79
914445	T914445	Circuit Breaker-2 amp	\$68.04
914480	T914480	Driven Roller Shaft	\$92.23

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
914482	T914482	Shaft	\$65.02
914491	T914491	Recirc. Pump-230V (M6B series)	\$1,431.33
914507	T914507	Bearing w/Groove (5000RA)	\$8.60
914568	T914568	Tube-Filter/Thermowell	\$32.76
914596	T914596	Filter EMI	\$326.84
914712	T914712	Dryer Front Cover (270RA)	\$124.17
914731	T914731	Orifice-0.25	\$50.90
914732	T914732	Orifice-0.125	\$9.42
914737	T914737	Chain-(62P;31 links)	\$121.40
914738	T914738	Chain (66P; 33 links)	\$176.72
914800	T914800	Quencharc	\$49.46
914813	T914813	Block, R.H.	\$1.95
914876*			
3H4056	T3H4056	Solution Heater (180LP)	\$720.20
914985	T914985	Block, L.H.	\$19.68
918211	T918211	Filter (460RA)	\$54.62
918641	T918641	Water Inlet Valve (M43)	\$84.48
918889	T918889	Orifice (.22 I.D.)	\$30.99
918925	T918925	Tube-35" x1.0" (ML300/180LP/270RA/3000RA)	\$87.32
918927	T918927	Tubing-35.75x.375	\$10.00
918928	T918928	Tubing-35.75x.375	\$33.20
918938	R918938	Developer Rack (270RA)	\$4,482.38
918943	R918943	Detector/Entrance X-Over Assy. (270RA)	\$928.79
918944	R918944	Dev/Fix X-Over Assy. (270RA)	\$1,823.46
943078*			
3H8738	T3H8738	Chain-Drive 32 links	\$122.03
943092	T943092	Tubing	\$9.28
943202	T943202	Pulley	\$31.74
943205	T943205	Tube, Replen.	\$27.78
943378	T943378	Varistor	\$9.91
943393	T943393	Chain, Dr. Motor	\$47.57
943402	T943402	Tube	\$13.99
943442	T943442	Filter, Solution	\$22.11
943457	T943457	5/8" Straight Tube (ML300/180LP/270RA/3000RA)	\$64.20
943468	T943468	Shaft, Drive (180 & 270RA series)	\$52.29
943471	T943471	Bushing (3000RA)	\$54.18
N.A.	TL943471	Bushing Assy. (3000RA)	\$54.18
943473	T943473	Shaft (180LP)	\$61.30
943629	T943629	Coupling, Body (Laser)	\$38.63
943665	T943665	Tubing	\$12.08
943666	T943666	Water Solenoid Valve	\$180.81
943667	T943667	Tubing, Formed	\$20.68
943700	T943700	Coupling (Male)	\$10.78
943701	T943701	Quick Disconnect	\$72.39
943712*			
550781	T550781	Gear, X-Over	\$10.03
953828	T953828	Spacer, Wash	\$5.00
955736*			
1C4593	T1C4593	Replenishment Pump (180,270,460,480,3000,5000RA)	\$1,808.15
955745*			

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
4E8210	T4E8210	Quad DC Power Supply	\$575.57
955595	T955595	Brkt-Stop, Non-Drive	\$2,041.46
963237	T963237	Spacer	\$9.28
963238	T963238	Spacer	\$10.00
963302	R963302	Developer Turnaround Assy. (M35A)	\$749.80
963541	T963541	Valve (1 gpm)	\$129.84
963545	T963545	Wash Solenoid Valve	\$227.56
963548	T963548	Dryer -Front Cover (180LP)	\$152.71
963611*			
9B8359	T9B8359	Blower Motor (M6B series)	\$990.53
963965	T963965	Idler, Roller	\$29.42
968064*			
3H9104	T3H9104	Support Casting-Dryer End	\$575.63
968067	T968067	Plastic Processing Tank (M6 seroes)	\$4,082.93
968199	T968199	Stud	\$72.07
968386	T968386	Main Drive Motor (M6I)	\$1,522.92
968389	T968389	AC Wire Harness Assy.	\$524.10
968433	T968433	Belt, Dryer (56.5")	\$188.37
968444	T968444	Spacer	\$102.19
968446	T968446	Spacer	\$20.68
968736	T968736	Shaft, Dryer (270RA & Integrated)	\$829.93
968781	T968781	Adapter-Harness	\$22.86
968802	T968802	800 PCB Board (M35 series)	\$622.89
968810*			
981299	T981299	Shaft, Dryer (180)	\$206.45
981345	T981345	Guard (M6 series)	\$67.18
981430	T981430	Roller	\$106.41
981719*			
TLUN5541	TTL5541	Manifold Test Kit-Vacuum (Laser)	\$173.65
981720*			
TL5542	TTL5542	Manifold Test Kit-Pressure (Laser)	\$3,247.00
981721*			
TL5544	TTL5544	Coupling-Red (Laser)	\$76.50
981722*			
TL5543	TTL5543	Coupling-Black (Laser)	\$88.00
981774*			
9B8359	T9B8359	Blower Motor (M6B series)	\$990.53
N.A.	T9870	Drip Pan w/ Drain (M6 series)	\$864.96
990010	T990010	Shoulder Screw	\$138.73
990019	T990019	Shaft, Idle (270RA wash)	\$295.85
990022	T990022	Shaft (180LP)	\$57.14
G9901918	T9901918	Plastic Key	\$7.11
990242	T990242	Guide Shoe, Large X-Over (smooth)	\$569.33
990243	T990243	Guide Shoe, Small X-Over (smooth)	\$110.50
990305	T990305	Guide Shoe, X-Over (M35 series)	\$100.11
990308	T990308	Guide Shoe, X-Over (M35 series)	\$434.51
990310	T990310	Spray Header Guide Shoe Assy. (M6I/460RA/480RA)	\$1,204.26
G9903246	T9903246	Pressure Gauge Vacuum (ML700)	\$48.38
990400	T990400	Fitting	\$9.15
990401	T990401	Quick Disconnect (Female)	\$75.79

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
990402	T990402	Recirc. Pump 230/240V	\$1,247.63
990415	T990415	Fitting	\$9.49
990416	T990416	Fitting, Recirc. Lrg.	\$13.10
990417	T990417	Nut, Tank, Small	\$10.72
990418	T990418	Nut, Gland	\$10.00
990419	T990419	Nut, Tank Bulk , Lrg.	\$10.00
990795	T990795	Isolation Relay	\$669.55
991743*			
1F5313	T1F5313	Solid State Relay	\$228.75
N.A.	Battery	Round Battery for DDTG1	\$30.63
N.A.	CHINABlack	X-Ray Film Marker (each) *Black*	\$3.46
N.A.	CHINARed	X-Ray Film Marker (each) *Red*	\$3.46
N.A.	CWK1	Cold Water Kit (M6AN/AW)	\$153.08
N.A.	CWK-PT	Cold Water Kit)for M6A w/T968067 tank)	\$205.42
N.A.	DF1	Drain Fitting	\$5.77
N.A.	DDTG1	Digital Dev. Temp. Gauge	\$42113
N.A.	DTA-M7	Developer Thermowell Assy.-M7	\$1,128.25
N.A.	DTA-M35110	Developer Thermowell Assy.-M35A/2000A (110-V)	\$1,400.00
N.A.	DTA-480RA	Developer Thermowell Assy.-480RA	\$1,388.49
N.A.	DTA-M6A	Developer Thermowell Assy.-M6A	\$3,586.26
N.A.	DTA-M6B	Developer Thermowell Assy.-M6B	\$1,722.32
N.A.	ELBORELIEF	Elbow Relief	\$3.07
N.A.	EZC1417	Clean-Up Film (14x17)	\$168.49
N.A.	FPK1	Front Plate Kit w/Thermowell and DDTG1	\$765.55
N.A.	TILFTPD	Illuminator Foot Peddle with 2-Way Switch	\$535.88
N.A.	IDProtector	I.D. Camera Protector	\$58.09
N.A.	M6ATAACK	M6A Tank Conversion Kit	\$898.76
N.A.	M6B.5amp	Circuit Breaker (M6B series-.5amp)	\$49.34
N.A.	M6DRAIN	Drain Valve Assy.	\$82.32
N.A.	T35MAMO	Dev. Rack for M35 Processor (Mammo processing)	\$2,713.90
N.A.	M6ADH-Assy.	M6 Dryer Heater Assembly	\$393.87
N.A.	M6BD/FDRAIN	Complete M6B D/F Drain Valve Assy.	\$140.99
N.A.	M6BMPH	M6B Main Power Harness	\$171.52
N.A.	MAMODEVKIT	Kit for Adjustable Shoes for M35A-M Dev. Rack Assy.(kit comes without adjustment tool)	\$89.96
N.A.	M7CHAIN	Replacement Chain for TI-BA M7B Main Drive Motor	\$30.62
N.A.	TM74WDPK	M7 4-wire Dryer Plug Kit	\$10.00
N.A.	PANSET3	Panel Set for M6 series Processors (Powder-coated - Greige)	\$2,600.71
N.A.	Pvcdrvalve	PVC 2-Way Ball Valve (1/2"x1/2")	\$31.74
N.A.	REMOVER	Utility Gear Remover	\$114.87
N.A.	RMAINDRIVE	Main Drive Assembly (M6AN/AW series)	\$1,516.63
N.A.	RMAINNEW	Main Drive Assy. (new M6B/460RA/480RA)	\$1,892.57
N.A.	RMAIN5000	Main Drive Assy. (5000RA)	\$2,043.97
N.A.	RMAIN270RA		\$ 493.30
N.A.	RMAIN35	Main Drive Assy. (M35)	\$616.90
N.A.	Speed-Up M6	60-Second Speed-Up Kit for M6	\$68.14
N.A.	Speed-Up 35	90-Second Speed-Up Kit for M35	\$98.79
N.A.	SH-8	Keypad Lock for Portable X-Ray	\$830.55
N.A.	TM1	Tube Mender (5/8")	\$6.13

<u>OEM Catalog #</u>	<u>TI-BA Catalog #</u>	<u>Description</u>	<u>List Price</u>
N.A.	TM6A-Therm	M6A Thermostat Assy. w/OT	\$3,128.68
NA.	TM35AM-L	<i>Non-Drive Side Shoe Spacer (Mammo Dev. Rack)</i>	\$10.00
N.A.	TM35AM-S	<i>Drive-Side Shoe Spacer (Mammo Dev. Rack)</i>	\$10.00
N.A.	TM35SST	Stainless Steel "T" Roller (M35 series)	\$277.13
N.A.	TPROTECTOR	Feed Tray Overlay Protector (M6 series)	\$61.22
N.A.	TM7PROTECT	Feed Tray Overlay Protector (M7 series)	\$61.22
N.A.	TR1	Tube Reducer (.75 to .5)	\$4.63
N.A.	UTILM35	Utility Roller (M35 series)	\$23.01

Multiloader Accessories

<u>OEM CAT #</u>	<u>TI-BA CAT #</u>	<u>Description</u>	<u>List Price</u>
N.A.	KUP3	Kover Up 3 (<i>Protective Cover for ML300</i>)	\$656.65
N.A.	KUP7	Kover Up 7 (<i>Protective Cover for ML700</i>)	\$802.92
N.A.	T700DECAL	Keyboard Face- <i>Decal Only</i>	\$1,123.48
TL1481	TTL1481	Alignment Tool	\$11.74
TL2170	TTL2170	Pinch Clamp	\$253.45
TL2431	TTL2431	Air Meter w/case	\$190.89
TL4391	TTL4391	Cable-Laptop Interface	\$158.70
TL4398	TTL4398	Small Chip Removal Tool	\$130.10
TL4488	TTL4488	Cassette Latch Tool	\$224.91
TL4722	TTL4722	Pressure Gauge (ML700)	\$47.31
TL4723*			
G9903246	T9903246	Pressure Gauge Vacuum (<i>ML700</i>)	\$65.65
TL4728*			
SP29010188	TSP29010188	Adapter Tool (<i>ML700</i>)	\$50.21
TL4729	TTL4729	Antistatic Wristlet (<i>ML700</i>)	\$46.91
TL4740*			
G9901918	T9901918	Plastic Key	\$7.11
TL4804	TTL4804	Gage-Cassette Opener	\$583.28
TL4830	TTL4830	Blow Pipe Positioner (<i>ML700</i>)	\$120.58
H4070061	T4070061	O-Ring (<i>ML300/270RA/3000RA/180LP</i>)	\$5.00
4070081	T4070081	O-Ring (<i>ML300/270RA/3000RA/180LP</i>)	\$ 6.86
4278842	T4278842	Screw (<i>M4x6</i>)	\$5.00
H4279240	T4279240	Screw (<i>M3x20</i>)	\$5.00
H4279346	T4279346	Round Pan Head Screw (<i>M5x16</i>)	\$5.00
4281551	T4281551	Screw (<i>M3x16</i>)	\$5.00
4281571	T4281571	Screw (<i>M6x12</i>)	\$5.00
4281756	T4281756	Screw (<i>M6x16</i>)	\$5.00
H4282691	T4282691	Screw	\$5.00
4282806	T4282806	Screw (<i>M4x16</i>)	\$6.96
4282836	T4282836	Screw	\$5.00
4283126	T4283126	Screw (<i>M4x10</i>)	\$5.00
H4480176	T4480176	Washer	\$5.00
4485810	T4485810	Washer	\$1.69
4500566	T4500566	Circlip	\$5.00
H4524181	T4524181	Relay, K1	\$31.94
4529176	T4529176	Plug	\$32.70
4532836	T4532836	Fuse Cap-60 Hz (<i>ML300</i>)	\$35.41
4533351*			
H4533356	T4533356	Microswitch	\$15.15
4533676*			
H4533671	T4533671	Fuse (US) 1.6A	\$10.00
H4534026	T4534026	Fuse (US) 4.0A	\$16.18
H4534036	T4534036	Fuse (US) 6.25A	\$10.00
4850311	T4850311	Nut	\$5.00
4853036	T4853036	Nut	\$5.00
4853111	T4853111	Nut (<i>M3</i>)	\$5.00
4853126*			
H4853121	T4853121	Nut (<i>M4</i>)	\$5.00

<u>OEM CAT #</u>	<u>TI-BA CAT #</u>	<u>Description</u>	<u>List Price</u>
4853256	T4853256	Nut	\$2.67
5001573	T5001573	Foam Tape (ML300/270RA/3000RA/180LP)	\$10.00
7094824	R7094824	Reconditioned Mutliloader Cassette (35x43 cm)	\$697.34
9152241	T9152241	Cap	\$118.63
9152271	T9152271	Cap	\$18.15
9180757	T9180757	Sensor B7/B8	\$87.76
H9181886	T9181886	Photocell	\$69.69
9182561	T9182561	Actuator	\$10.44
H9183486	T9183486	Inductive Sensor	\$538.97
9184847*			
H9304840	T9304840	Magazine Sucker (4/pkg.)	\$183.27
H9184396	T9184396	Ring, Foam Rubber (XML300/700)	\$10.00
H9185636	T9185636	Timing Belt (ML700)	\$46.31
H9185646	T9185646	Timing Belt	\$51.98
H9185676	T9185676	Cassette Sucker (ML300/700)	\$75.22
9186461	T9186461	Fan-60 Hz-Frame (ML700)	\$278.40
9187246*			
H9187241	T9187241	Air Filter	\$22.04
9187631*			
H9182216	T9182216	Photo Sensor Assy. (ML700)	\$611.11
9187851*			
H9187856	T9187856	Air Filter	\$12.08
H9187903	T9187903	Fan-(50-60 Hz) (ML700)	\$347.89
9188568*			
H9188586	T9188586	Belt-Tunnel (ML700)	\$32.82
H9188576	T9188576	Round Belt (ML700)	\$40.26
9189907	T9189907	Infill Panel, RH	\$337.87
9189947	T9189947	Infill Panel, LH	\$427.20
9189996	T9189996	Label, Film Size	\$39.88
H9190136	T9190136	Ribbon Cable (ML700)	\$163.74
H9190306	T9190306	Ribbon Cable Assy. (ML700)	\$246.96
H9192081	T9192081	Wire Assy.	\$42.40
H9192226	T9192226	Cable Proc Start (ML700)	\$186.54
H9193661	T9193661	Foam Rubber (12x12 10M)	\$54.56
H9194556	T9194556	Reflective Foil (180 labels)	\$55.50
N.A.	S9194561	Reflective Label for T9284356	\$80.20
9198033*			
9198036	T9198036	Water Supply Assy. (ML300/ML700)	\$178.50
9198236	T9198236	Humidifier Fan (ML300/ML700)	\$267.18
H9199286	T9199286	Magnetic Valve	\$228.31
9211891	T9211891	Spacer	\$5.00
9212941	T9212941	O-Ring-Film Pocket	\$5.00
9218836	T9218836	Ribbon Cable (ML300)	\$172.87
9220210	T9220210	Screw	\$16.58
9221160	T9221160	Belt	\$49.01
9221310	T9221310	Belt (ML300)	\$35.03
9221400	T9221400	Roller Shaft for Daylight System (ML300)	\$163.30
H9223151	T9223151	Belt-Timing (ML300)	\$21.39
H9223400	T9223400	Timing Belt (ML300)	\$19.11
H9223410	T9223410	Spring (XML300/ML700)	\$2.56

<u>OEM CAT #</u>	<u>TI-BA CAT #</u>	<u>Description</u>	<u>List Price</u>
9224840*			
H9304840	T9304840	Magazine Sucker (4/pkg.)	\$183.27
H9226140	T9226140	Belt	\$32.63
H9226490	T9226490	Manometric Switch	\$444.72
H9227030	T9227030	Bearing	\$10.00
9227320	T9227320	Label	\$33.08
H9228690	T9228690	Relay, K1	\$328.99
H9228920	T9228920	Fan (ML300)	\$342.97
H9228980	T9228980	Sensor B24, B60	\$545.96
H9228991	T9228991	Sensor B20 (ML700)	\$75.98
H9229880	T9229880	Spring (XML300/ML700)	\$5.00
H9281906	T9281906	Power Supply (ML700/700+)	\$1,993.88
H9284026	T9284026	Sensor-B14 (ML700)	\$153.66
H9284356	T9284356	B6 Sensor Kit (ML700)	\$761.20
H9285601	R9285601	Rebuilt ML700 Vacuum Pump (60 Hz)	\$1,296.72
H9304840	T9304840	Magazine Sucker (4/pkg.)	\$183.27
N.A.	BracketML	Bracket for T9284356 B6 Sensor Kit	\$40.88
N.A.	SENSOR700	Sensor Bracket (for Kit) (ML700)	\$87.63
N.A.	TRANASSY	QC Transistor Assy. for T9284356 B6 Sensor Kit	\$146.15

Dry Laser Accessories

<u>OEM CAT #</u>	<u>Description</u>	<u>List Price</u>
SP1C8261	Isopropyl Alcohol Cleaner	\$32.01
1C8287	Alcohol Wipes (50/pk) Lint Free	\$85.18
TL5303	Gloves Neophene-Acid Resist	\$15.56
TL5568	7.5M PT-PT Network Cable	\$ 139.47
TL5574	Digital Temp. Meter Kit 11 Probe	\$ 708.34
162786*		
<i>SPK162786</i>	E-Ring, Retaining	\$5.00
3E5676*		
<i>34704189397</i>	Paper Core Drum Stand	\$.61
SP22000111249	Static Bag	\$2.27
SP26101211303	Oil, Krytox GPL 106, 1/2-oz. syringe	\$ 104.98
SP26101167885	Filter-Air, Right Side	\$83.78
SP34740189397*		
<i>34704189397</i>	Paper Core Drum Stand	\$.61
44002271821	8100 PM Kit	483.11
SP6E7324	Felt Insulator-Mount Cast. Vers.	\$6.00
SP66000001546*		
<i>SP8E9166</i>	Kit-Stripper Replacement Roller Bearing	\$ 38.53
SP74040183209	Stripper Spring Replacement Kit	\$ 3.27
74040163318*		
<i>SP74040183209</i>	Stripper Spring Replacement Kit	\$3.27
74040179728*		
<i>SP74040183209</i>	Stripper Spring Replacement Kit	\$3.27
SP74040178944	Gasket Kit	\$80.83
SP74040184439	Felt Replacement Kit	\$27.08
SP74040189669	8100 Filter	\$ 355.26
SP74999991255	Felt Plate Assy.	\$ 56.69
SP78801828027	Webril Bag of Handi-Pds 4x4 (100% cotton)	\$18.48
78803374756*		
<i>1C8261</i>	Isopropyl Alcohol Cleaner	\$ 21.51
SP78805340599	Cable Video External (10M)	\$162.21
SP78805346950	Cable Video External	\$ 107.22
SP78805346976	Cable Image External	\$ 323.07
SP78807526054	Cable Assy-Fiber Opt.	\$205.62
SP78807740978	3-Meter Keypad Cable Assembly 969	\$481.14
SP78805341340	Cable Image Ext 10 Mtr. Shielded	\$ 460.53
SP78807741356	UKEIB Assy (HQ 969 Univ. KP Adapt. Box)	\$4,615.40
SP78807742016	DEIB Ay. Not Siemens Modality	\$3,701.98
SP78807907478	EVEIB Assy <i>w/exit Sync.</i>	\$9,068.57
78809240035*		
<i>SP8E2644</i>	Insulation-TOP	\$396.04
SP78809451293	990 Switch Assy. (Exit)	\$106.94
SP78809451301*		
<i>SP8E2330</i>	Assembly-Switch, Optical Flag	\$ 35.99
78809456508*		
<i>SP6E7324</i>	Felt Insulator-Mount Cast. Vers.	\$6.00

<u>OEM CAT #</u>	<u>Description</u>	<u>List Price</u>
SP78809491323	(693CO V5.0) 32MP x 9 Bit Input Mod	\$13,207.92
SP78809459791*		
<i>SP8E9280</i>	Filter-Air, Right Side Panel	\$20.95
SP8E8545	AY-Power Supply	\$2,167.72
SP78809490549*		
<i>SP8E7331</i>	Charcoal Filter (8700)	\$ 258.19
SP78809594951	Filter, Air, Power Module	\$ 4.46
78800537793	Cloth (150/pk) Lint Free	\$250.09
SP78809594951	Kit-Stripper Replacement Roller Bearing	\$ 53.37
SP78809632041	Kit-8300 Charcoal Filter	\$331.51
SP78809632199	Upper Felt Pad AY (for stripper)	\$123.77
SP78810005229	Kit-8300 Faz Filter	\$ 103.09
78811320692	Lower Felt Pad Assy (for stripper)	\$ 124.41
78811321047	8300 Idler Roller	\$ 303.88
SP78865640011	Retaining Ring	\$.27
SP78999828680	One Bag-with Twist Tie	\$1.19
SP78999829159	10,000 Cycle PM Kit (8700)	\$353.71
SP8E2330	Assembly-Switch, Optical Flag	\$ 35.99
SP8E2402	Kit-Optics Module, Gen 3 (exchange item)	\$19,342.24
SP8E2643	Insulation-BOTTOM	\$414.38
SP8E2644	Insulation-TOP	\$396.04
SP8E7286	Drum Assy-8700 (exchange item)	\$4,604.60
SP8E7331	Charcoal Filter (8700)	\$ 258.19
SP8E9166	Kit-Stripper Replacement Roller Bearing	\$ 38.53
SP8E9280	Filter-Air, Right Side Panel	\$20.95
SP96000000397	Barrier-Air	\$106.33
96000000439*		
<i>SP8E2643</i>	Insulation-BOTTOM	\$ 414.38
SP96000000578	O-Cello Pad	\$ 4.25
SP96000000644	Stainless Steel Clean & Polish	\$ 21.36
SP96000000669	Troubleshooter Cleaner	\$25.06
SP96000001188	Tack Cloth	\$ 4.30
SP96000003432	Roller Assy.	\$ 88.84
SP96000004042	Pressure Roller Rulon 0406 Flat	\$ 4.54
SP96000016590	Oil-Dow #710, Silicone 2 (2-oz.)	\$46.21
SP96000027555	Filter-Door	\$21.54

RACK KITS

Rack Parts Kits **contain the most-frequently replaced parts used to rebuild the rack.**

**** SIDEPLATES AND GUIDE SHOES ARE NOT INCLUDED ****

<u>TI-BA</u> <u>Catalog Number</u>	<u>Description</u>	<u>List Price</u>
TD35A	M35 <i>Dev.</i> Rack Parts Kit	\$ 1,685.52
TF35A	M35 <i>Fixer</i> Rack Parts Kit	\$ 1,809.43
TW35A	M35 <i>Wash</i> Rack Parts Kit	\$ 2,033.59
TDRK7	M7 <i>Dev.</i> Rack Parts Kit	\$ 1,939.86
TFRK7	M7 <i>Fixer</i> Rack Parts Kit	\$ 1,718.32
TWRK7	M7 <i>Wash</i> Rack Parts Kit	\$ 2,005.96
TDRK6A	M6 <i>Dev.</i> Rack Parts Kit	\$ 2,094.03
TFRK6A	M6 <i>Fixer</i> Rack Parts Kit	\$ 1,731.56
TWRK6A	M6 <i>Wash</i> Rack Parts Kit	\$ 1,517.49
TWRK6/AW	M6AW <i>Wash</i> Rack Parts Kit	\$ 1,947.77

SUPER RACK KITS

Super Rack Kits contains *all replacement parts* as contained in **Rack Parts Kit** above.
SIDEPLATES ARE INCLUDED - ** GUIDE SHOES ARE NOT INCLUDED **
 (nor does kit contains everything needed to build
 a complete rack)

TI-BA		
<u>Catalog Number</u>	<u>Description</u>	<u>List Price</u>
TD35SRK	M35 Super Developer Rack Kit (includes sideplates)	\$ 1,895.41
TF35SRK	M35 Super Fixer Rack Kit (includes sideplates)	\$ 1,937.29
TW35SRK	M35 Super Wash Rack Kit (includes sideplates)	\$ 2,246.21
TDM7SRK	M7B Super Developer Rack Kit (includes sideplates)	\$ 2,109.26
TFM7SRK	M7B Super Fixer Rack Kit (includes sideplates)	\$ 1,845.65
TWM7SRK	M7B Super Wash Rack Kit (includes sideplates)	\$ 2,177.30
TD6ASRK	M6A Super Developer Rack Kit (includes sideplates)	\$ 2,324.02
TF6ASRK	M6A Super Fixer Rack Kit (includes sideplates)	\$ 1,928.87
TW6ASRK	M6A Super Wash Rack Kit (includes sideplates)	\$ 1,731.53
TW6AWSRK	M6AW Super Wash Rack Kit (includes sideplates)	\$ 2,160.37

Transport Systems

TI-BA Catalog Number	Description	List Price
Transport-35	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash</i>)	\$9,270.35
Transport-35M	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash</i>)	\$9,409.43
Transport-MinR	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash</i>)	\$9,893.41
Transport-M7B	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash</i>)	\$13,778.89
Transport-AN	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash/Squeegee</i>)	\$11,410.90
Transport-AW	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash/Squeegee</i>)	\$12,489.04
Transport-M6B	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash/Squeegee</i>)	\$12,998.96
Transport-270RA	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash & Dryer</i>)	\$17,236.06
Transport-3000RA	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash & Dryer</i>)	\$15,561.99
Transport-480RA	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash</i>)	\$13,369.44
Transport-5000RA	Transport System (<i>consists of Entrance,/2 X-Over/Dev/Fix/Wash</i>)	\$12,798.56

TUBE KITS

Tube Kits contain all parts needed to *re-tube* the referenced processor.




<u>TI-BA Catalog Number</u>	<u>Description</u>	<u>List Price</u>
Tube- M35A	Tube Kit for M35A Processor	\$ 324.99
Tube- M7B	Tube Kit for M7B Processor	\$ 499.65
Tube- AN	Tube Kit for M6AN Processor	\$ 416.37
Tube- AW	Tube Kit for M6AW Processor	\$ 449.67
Tube- M6B	Tube Kit for M6B Processor (<i>S/N 9,999 and below</i>)	\$ 516.37
Tube- M6B1	Tube Kit for M6B Processor (<i>S/N 10,000 and above</i>)	\$ 770.59
Tube- 460	Tube Kit for 460RA Processor	\$ 699.48
Tube- 180	Tube Kit for 180LP/LPS Processor	\$ 782.77
Tube- 270	Tube Kit for 270RA Processor	\$ 649.57
Tube- 480	Tube Kit for 480RA Processor	\$ 449.67
Tube- 3000	Tube Kit for 3000RA Processor	\$ 584.38
Tube- 5000	Tube Kit for 5000RA Processor	\$ 866.05

Tank Fitting KITS


Kits contain all small parts on respective tanks. Tanks, weirs, heat exchangers, tank clips, tank seals sold separately.

<u>TI-BA Catalog Number</u>	<u>Description</u>	<u>List Price</u>
TFix-Kit	Tank Fitting Kit (<i>for M35 and M7B Fixer Tanks</i>)	\$146.54
TWash-Kit	Tank Fitting Kit (<i>for M35 and M7B Wash Tanks</i>)	\$ 91.58
TDev35-Kit	Tank Fitting Kit (<i>for M35 Developer Tanks</i>)	\$ 91.58
TDev7B-Kit	Tank Fitting Kit (<i>for M7B Developer Tanks</i>)	\$ 56.02

TRANSPORT SYSTEM CHART

DESCRIPTION	M20	M35A	M35A-M	2000/A	MIN-R	M7B	180LP
Detector Entrance	R8B66682	R5B6163	R5B6163	R5B6163	T3E1945	R474481	R8B6682
Developer	R474483	R261275	R7C7862 Or T35Mamo	5E2458	R3E1944	R474483 S474483	R7C8328  (Phase 4)
Dev/Fix Crossover	R335975	R333795	R333795	R333795	T3E1943	R335975 (Rubber) Or S335975 (Silicone)	R918944 R8B6576 (S/N: 3,000-13,000)
Fixer	R495748	R748165	R748165	R748165	R748165	R495748	R8B6679
Fix/Wash Crossover	R335976	R261276	R261276	R261276	R261276	R335976	R9B8445 R9B8832 (S/N: 3,000-13,000)
Wash	R240682 (Silicone)	R261277	R261277	R261277	R261277	R495749	R8B6612 R5B6575 (S/N: 3,000-13,000)
Squeegee	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.
Turnaround Assy. (F/W)	N.A.	R748164	R748164	N.A.	3E1946	R333789	N.A.
Turnaround Assy. (D/F)	N.A.	R963302	914265	N.A.	N.A.	R335977	
Dryer Rack	N.A.	333796	333796	333796	N.A.	R8B7075	R9B8086 
Exit Rack	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.	R8B6577 

NOTE:  MUST USE DETECTOR ASSEMBLY #R8B6682 WITH THIS PART # DEV. RACK.

 Requires a Trade-In or Pay TI-BA's selling price!

DEV/FIX X-OVER PART # R918944 (USE TROUGH #914920)

DEV/FIX X-OVER PART # R8B6576 (USE TROUGH #1C7229)

FIX/WASH X-OVER PART #9B8445 (USE TROUGH #914922)

FIX/WASH X-OVER: PART #R9B8832 (USE TROUGH #1C7230)

180LPS (WITH SERIAL NUMBERS ABOVE 13,000) USE SAME RACK CONFIGURATION AS 180LP (WITH SERIAL NUMBERS ABOVE 3,000)

TRANSPORT SYSTEM CHART

DESCRIPTION	270RA "Stand Alone"	270RA "INTEGRATED"	3000RA	M6AN	M6AW	M6B	460/480RA	5000RA
Detector Entrance	R8B6682	R8B6682	R8B6682	R1C7377 (3-roller) R240683 (2-roller)	R1C7377 (3 roller) R240683 (2 roller)	R1C7377 (3 roller) R240683 (2 roller)	R1C7377 (3 roller) R240683 (2 roller)	R1C7377 (3 roller) R240683 (2 roller)
Developer	R7C8328 ♦ (Phase 4)	R7C8328 ♦ (Phase 4)	R7C8328 ♦ (Phase 4)	R240680	R240680 S240680	R240680	R240680	R8B7293
Dev/Fix Crossover	R918944	R918944	R5B6576	R474478 (Rubber)	R474478 (Rubber)	R240684 (Silicone) S240684 (a mammo option)	R240684	R8B6965
Fixer	R8B6679	R8B6679	R8B6679	R240681	R240681	R240681	R335534	R8B7294
Fix/Wash Crossover	R9B8445	R9B8445	R9B8832	R779155	R779155	R779155	R261897	R3H8726
Wash Rack Assy.	R8B6612	R8B6612	R5B6575	R512893 (Rubber)	R240682 (Silicone)	R240682 (Silicone)	R240682	R8B6969
Squeegee	N.A.	N.A.	N.A.	R240686	R240686	R240686	R240686	R8B6970
Turnaround Assy. (Phenolic "B" Roller)	N.A.	N.A.	N.A.	R333789	R333789	R333789	N.A.	N.A.
Turnaround Assy. (Rubber "B" Roller)	N.A.	N.A.	N.A.	R335977	R335977	R335977	N.A.	N.A.
Dryer Rack	R696093	R9B8436	R5B6580	N.A.	N.A.	N.A.	N.A.	N.A.
Exit Rack	R696092	N.A.	R5B6579 ★	N.A.	N.A.	N.A.	N.A.	N.A.

NOTE: ♦ MUST USE DETECTOR ASSEMBLY #R8B6682 WITH THIS PART # DEV. RACK.

DEV/FIX X-OVER PART # R918944 (USE TROUGH #914920)

FIX/WASH X-OVER PART #9B8445 (USE TROUGH #914922)

180LPS (WITH SERIAL NUMBERS ABOVE 13,000, USE SAME RACK CONFIGURATION AS 180LP (WITH SERIAL NUMBERS ABOVE 3,000)

★ Requires a Trade-In or Pay TI-BA's Selling Price!

DEV/FIX X-OVER PART # R8B6576 (USE TROUGH #1C7229)

FIX/WASH X-OVER: PART #R9B8832 (USE TROUGH #1C7230)

Pumps and Motors


Machine Model	Replenishment Pump	Recirculation Pump	Main Drive Motor	Blower Motor Assy.
M8	T594394	T594391		
M20	T651861 (230V)	T418702 (230V)	482241 (230V) **Discontinued**	482244 (230V)
M35A/AM & 2000A	T240227 (115V)	990403 (115V)	779519 (M35A-M) 246839 (110V) (M35A & 2000)	T240225 (115V)
M35/m & 2000	T651861 (230V)	T333519 (230V)	651863 (230V)	651864 (230V)
Min-R	T651861 (230V)	T990402 (230V)	651863 (230V)	T2E4249 (230V)
M7/A/B	T594760 (110V)	T418709 (110V)	8B6920 (110V)	R532237 (110V)
180LP	T1C4593 (230V)	1C4594 (230V)	943105	748638
270RA	T1C4593 (120V)	1C4594 (230V)	779325 (Gear Box) 914547 (Motor)	748638
3000RA	T1C4953 (120V)	1C4594 (230V)	1C4591 (Gear Box) 1C4589 (Motor)	9B8251
5000RA	5B6353	T5B6352	5B6376	8B6625
M6AN/AW	T914429 (230V)	T418702 (230V)	T3H4475 (60 Hz.)	T748683 S748683 (50 Hz)- European Version T474694 (w/pulleys)
M6B (serial numbers below 10,000)	T914429 (230V)	T622017 (230V)	T655658 (230V)	T748683 (230V)
M6B (serial numbers above 10,000)	T333659 (230V)	T914491 (230V)	T333658 (230V)	T9B8359 (230V) S9B8359 (50 Hz)- European Version
460RA	943211 (230V)	T914491 (230V)	779325 (Gear Box)- 460RA	911428 (230V)
480RA			914547 (Motor)-460RA 779325 (Gear Box)- 480RA 955703 (Motor)-480RA	T9B8359 (230V)

Didn't Find a Particular Item or Part Number?

Don't Despair!

We have an un-ending supply of used equipment & parts . . .
If you don't see a part or item listed, give us a call!

It may have been discontinued by the
OEM, but we often can find it here in a
"used" state . . .
. . . which may help get you out of a bind . . .

. . . could make you shine like a  . . .
. . . a hero . . .
in your customer's eyes . . .

It never hurts to ask . . .

. . . and it always **HURTS LESS!** when the invoice arrives . . .

Prices effective 10.15.10

Part number	Description	List Price	ACL2i	ACL2	ACL4	PoC 120/140	PoC 260	CR360
AC000110*								
8002396	CLNG SCREEN TRAY FOR POC <i>(with autoloop)</i>	60.00						
AC000111*								
8002412	ROLLER CLNG SCREEN POC <i>(with autoloop)</i> <i>(to be used with Blue Tray)</i>	100.00				X	X	
AC085003*								
AS000425	PoC Packing Crate	534.00			X	X	X	
AC097601*								
8002420	CLEANING SCREEN TRAY POC <i>w/o autoloop</i>	250.00						
AC097602*								
8002438	ROLLER CLNG SCREEN POC <i>w/o autoloop</i> <i>(to be used with Black Tray)</i>	100.00		X	X	X	X	
AK000201	GOP LICENSE CODE	1,040.00		X	X	X	X	
PR000018*				X	X	X	X	
PR000290	LABEL FOR CASS. 14"X17" POC AGFA	25.00		X	X	X	X	
PR000020	LABEL FOR CASS. 10X12" POC AGFA	12.00		X	X	X	X	
SK000002	AUTOLOOP KEY AY POC 140	1,365.00				X		
SK000003	SK AUTOLOOP UPG POC 140 SW 2.1.2	3,108.00				X		
SK000007	AUTOLOOP KEY AY POC 260	2,246.67					X	
SK000008	AUTOLOOP UPG POC260 SW 2.1.2	4,938.00					X	
SK000010	MCD POC260 AUTOLOOP-EXCH	2,664.00					X	
SK000012	SW V2.4.2 POC 260	8.00					X	
SK000013	VER 2.4 S/W FOR CR260	8.00					X	
SK000014	AUTOLOOP UPG POC 140 SW V2.4.2	2,877.00				X		
SK000015	AUTOLOOP UPG POC 260 2.4.2	4,486.68					X	
SK000019	OPTIC HEAD ACL2i	6,409.94	X					
SK000035	SP 3A+BOOT V57	45.62			X			
SK000036	SENSOR-CASSETTE PRESENCE	73.02		X	X	X	X	
SK000037	AUTOLOOP UPG POC 140 SW 2.5	6,846.42				X		
SK000038	AUTOLOOP UPG POC 260 SW 2.5	6,846.42					X	
SK000046*								
SK000096	CONTROLLER BOARD-EXCH	1,800.00						X
SK000047	IMCS BOARDS ASSY-EXCH	4,595.00						X
SK000048	3.SK LASER BOARD CR 360	245.88						X
SK000049	CCFL INVERTER ASSY-EXCH	2,436.00						X
SK000050	SK CCFL DRIVER BOARD CR 360	764.09						X
SK000051	ERASING LAMPS ASSY-EXCH	3,260.00						X
SK000052*								
SK000088	POWER SUPPLY CR360-EXCH	2,898.12						X
SK000053	OPTICAL HEAD ASSY-EXCH	7,700.74						X
SK000054	PM HEAD ASSY CR 360-EXCH	8,315.54						X
SK000055*								X
SPSK000321	HAYDON LEAD SCREW ASSY.	2,934.00						X
SK000056*								
SK000087	SK BARCODE CR 360	685.20						X
SK000057	SK 24V FILTER BOARD CR 360	198.24						X
SK000058	SK ENCODER BOARD CR 360	243.96						X
SK000059	SK FLEX PRINT CONNECTOR CR 360	196.80						X
SK000060	SK ROLLERS REPLACE KIT CR 360	1,569.64						X
SK000061	SK ROLLER MOTOR COUPLING CR 360	171.00						X
SK000062	SK LINEAR MOTOR COUPLING CR 360	222.48						X
SK000063	SK PLATE CARRIAGE ASSY CR 360	620.68						X
SK000064	SK ROLLERS MOTOR ASSY CR 360	1,575.42						X

SK000065	SK LOADER ASSY CR360	820.00							X
SK000066	SK LINEAR MOTOR CR 360	472.00							X
SK000067	SK TRAY ASSY CR 360	3,340.60							X
SK000068	SK ROOF ASSY CR 360	1,500.00							X
SK000069	SK POWER SW & LINE FILTER CR 360	303.78							X
SK000070	SK FUSES CR 360	606.24							X
SK000071	SK SENSOR J502 CR 360	179.28							X
SK000073	SK SENSOR J504 CR 360	179.28							X
SK000074	SK SENSOR J505 CR 360	174.78							X
SK000075	SK SENSOR J508 CR 360	179.08							X
SK000076	SK SENSOR J510 CR 360	147.36							X
SK000077	SK SENSOR J512 CR 360	151.86							X
SK000078	SK SENSOR J514 CR 360	179.28							X
SK000079	SK SENSOR J515 CR 360	191.40							X
SK000080	SK SENSOR J516 CR 360	201.30							X
SK000081	SK TOP COVER ASSY CR 360	1,600.54							X
SK000082	SK BOTTOM COVER ASSY CR 360	1,477.68							X
SK000083	SK ERASING LAMPS THERMO CR 360	154.68							X
SK000084	SK FUSES CR 360	209.64							X
SK000085	PM BOARD CR 360-EXCH	3,248.86							X
SK000086	SK SENSOR J513 CR 360	163.86							X
SK000088	SK POWER SUPPLY (LN FAN) CR 360	2,898.12							X
SK000099	SK SCAN HEAD ALIGN GAUGE CR 360	372.80							X
SK000100	TECH PLUG CR 360	203.28							X
SK000105	GOP DONGLE PARRALLEL + Lic-exch	2,808.00			X		X	X	
SK000106	GOP DONGLE PARALLEL + LIC-Exch	2,808.00			X		X	X	
SK000107	SOLENOID FOR AUTOLOOP CR140	590.69					X		
SK000108	AUTO-LOOP SOLENOID	637.50						X	
SK000114	W0 SENSOR IR FILTER POC 140	257.22					X	X	
SK000116	LOADER ASSY POC260-EXCH	2,433.60						X	
SK000117	LOADER ASSY ACL 4 / POC140	1,756.40			X		X		
SK000118	LOADER ASSY POC140	2,515.50					X		
SK000129	LINEAR MOTOR POC260-EXCH	1,602.00						X	
SK000133	SK FRONT LABEL POC 140 IDEXX	54.73					X		
SK000146	SK 4 RED ROLLERS PROC 120/140	1,508.98					X		
SK000155	ROLLER SPRINGS 50MM(POC) (set of 2)	90.12					X	X	
SK000159	INVERTERS MLX CON (POC120/140)	2,436.00					X		
SK000160	ERASE LAMP MLX CON (POC120/140)	2,436.00					X		
SK090278	L-BRACKET-UPPER	25.00	X	X	X		X	X	
SK200000	USB BOARD	2,913.00			X				
SK200001	KIT FOR MOTION MED-EXCH	1,665.60			X				
SK200002	SERVICE KIT ENHANCED MOTION MED	477.54			X				
SK200003	SENSOR BOARD	229.86			X				
SK200006	LASER CARD	315.16			X				
SK200007	5B/12V POWER SUPPLY	289.80			X				
SK200009	12V POWER SUPPLY (ROT)	580.00			X				
SK200012	STEP MOTOR SLIDE	327.24			X				
SK200014	ROLLER MOTOR MED KIT	132.00							
SK200016	INVERTER	449.58			X				
SK200019	SERIAL EPROM SCV/KIT	158.40			X				
SK200020	SENSOR PLUS WRIING	169.68			X				
SK200021	ROLLER MOTOR	984.18			X				
SK200023	12V POWER SUPPLY (ROT)	451.73			X				
SK200025	PLATE GUIDE SPRING	118.74			X				
SK200026	MOTION BOARD	1,662.92			X				

SK200027	PLATE GUIDE	398.52		X	X		X
SK200028	USB BOARD, ACL2I-EXCH	3,226.50		X	X		
SK200029	MOTION BOARD ACL2i	1,645.50	X				
SK240000	USB BOARD-EXCH	3,145.50			X		
SK240001	MOTION BOARD-EXCH	1,111.00			X		
SK240002	ERASE LAMPS BOARD	379.00			X		
SK240005	SENSOR BOARD	236.34			X		
SK240006	PM PLUS PM BOARD-EXCH	5,659.25			X		
SK240007	LASER BOARD	126.00			X		
SK240008	LINEAR WAYAY-ACL4/POC=EXCH	5,677.00	X		X	X	
SK240013	STEP MOTOR SLIDE	1,221.76			X		
SK240019	15V POWER SUPPLY	784.00			X		
SK240020	SENSOR PLUS WIRING	119.28			X		
SK240021	ERASE LAMPS BOARD	911.31			X		
SK250001	MOTION BOARD	1,631.88			X		
SK250002	SENSOR BOARD	226.38			X	X	X
SK250004	LASER BOARD	167.94			X	X	
SK250005	ROLLER MOTOR	786.29			X	X	X
SK250006	LINEAR MOTOR	354.18			X	X	
SK250007	LOADER MOTOR	348.54			X	X	X
SK250008	POWER SUPPLY AY-EXCH	2,971.00			X	X	
SK250009*							
SK250010	LOCATING SENSOR J502	198.00			X		
SK250010	LOCATING SENSOR J502	198.00			X	X	
SK250011	LOCATION SENSOR J505/J506/J508	221.88			X	X	
SK250012	LOCATING SENSOR J507	70.98			X		
SK250014	LOCATING SENSOR J513/J514	156.00					X
SK250015	R-ORIGIN SENSOR J515	126.00			X	X	X
SK250016	L-ORIGIN SENSOR J516	67.62			X	X	X
SK250017*							
SK25003	ERASE LAMPS SENSOR	524.58			X		
SK250018	ERASE LAMPS-EXCH	2,436.00			X		
SK250019	INVERTER-EXCH	1,814.00			X		
SK250020	SERIAL EPROM	229.86			X		
SK250021	COVER	2,459.82	X	X	X		
SK250022	USB BOARD-EXCH	3,184.50			X		
SK250023	PM ASSEMBLY-EXCH	5,710.75			X	X	
SK250024	OPTICAL HEAD-EXCH	5,543.25			X	X	
SK250025	ROTATION MOTOR DRIVE_EXCH	624.00			X	X	
SK250026	FUSES SET	64.56			X	X	
SK250027	FLEX CABLE 1300MM	82.86			X	X	
SK250028	FLEX CABLE LED PANEL	60.00			X	X	X
SK250029	FLAT CABLE MOTION TO USB	102.60			X		
SK250030	FLAT CABLE SENSOR TO MOTION	25.00			X	X	X
SK250031	CABLE MOTION TO INVERTERS	116.64			X	X	
SK250032	OPTICAL HEAD ST	5,409.18			X		
SK250033	ROLLER MOTOR COUPLING	180.00			X	X	X
SK250034	STEPPER MOTOR COUPLING	264.23			X	X	
SK250035	MAIN GROUND CONNECTION CABLE	298.68			X	X	X
SK250036	SCREW SET	32.83			X	X	X
SK250037	5V/12V POWER SUPPLY	355.80			X	X	X
SK250038	15V POWER SUPPLY	837.71			X	X	X
SK250039	12V POWER SUPPLY	820.07			X	X	
SK250040	CABLE ERASE LAMP SENSORS TO INVERT	25.00			X	X	
SK250041	CABLE LASER TO PM	25.00			X	X	

SK250042	COVER-BACK,METAL	155.04			X	X	X
SK250043	COVER-BOTTOM,METAL	148.50			X	X	X
SK250044*							
SK000155	ROLLER SPRINGS 50MM (POC) (pkg of 2)	90.12					X
SK250045	TIMING BELT	129.30				X	X
SK250046	ROLLERS KIT	1,875.48				X	X
SK250048	BOOT VER 4.2.1.53	142.42			X		
SK250049	BOOT VER 4.2.1.55	114.00			X	X	
SK250051	SENSOR J513 ASSY	429.97			X	X	
SK250052	SENSOR J514 ASSY	183.90				X	X
SK250053	BOOT VER 4.2.1.54	142.42			X		
SK250054	BOOT VER 4.2.1.56	112.50			X		
SK250058	FLEXIBLE PLATE GUIDE	172.81			X	X	X
SK255001	USB BOARD POC120/140-EXCH	3,440.00				X	
SK255003	ERASE LAMPS SENSOR	524.58				X	X
SK255004	COVER	2,635.98				X	X
SK255005	FLAT CABLE MOTION TO USB	45.50				X	
SK255006	MOTION POC 120/140 AUTO LOOP	1,008.00				X	
SK255018	ERASE LAMP POC120/140-EXCH	2,436.00				X	
SK260001	USB BOARD POC260 – EXCH	3,241.50					X
SK260003	X ENCODER BOARD	275.73					X
SK260004	MOTOR FLEX PRINT	209.70					X
SK260005	LASER BOARD	656.52					X
SK260006	PM ASSEMBLY – EXCH	7,056.25					X
SK260007	OPTICAL HEAD – EXCH	7,607.00					X
SK260008	LINEAR MOTOR	1,427.91					X
SK260009*							
SK000129	LINEAR BL MOTOR COUPLING	1,602.00					X
SK260010	ROTATION BL MOTOR DR-EXCH	2,619.00					X
SK260011	24V POWER SUPPLY	838.64					X
SK260012	LOCATING SENSOR J509	174.54					X
SK260013	LOCATION SENSOR J505/ J508	270.54					X
SK260014*							
SK000171	ERASE LAMPS SENSOR (PoC260)	303.08					X
SK260015	ERASE LAMPS-EXCH	3,140.00					X
SK260016	INVERTER ASSY-EXCH	2,436.00					X
SK260018	POWER SWITCH LINE FILTER	284.76					X
SK260019	FUSES SET	102.42					X
SK260020	BAR CODE READER-EXCH	2,762.64					X
SK260021	FLAT CABLE MOTION TO USB	37.83					X
SK260022	CABLE MOTION TO ERASE LAMP BOARD	39.39					X
SK260023	CABLE ERASE LAMPS SENSORS TO INVERT	37.83					X
SK260024	CABLE USB-PMT-LASER	337.46					X
SK260025	LINEAR ASSY POC260-EXCH	7,131.00					X
SK260026	TRAY AY-EXCH	2,710.50				X	X
TLUN8454977	SCREEN CLEANER	21.44	X	X	X	X	X
UNSP869623	LABEL- MOD	30.00	X	X	X	X	X
ST095870-02*							
PR000021	BAR CODE LABEL CASS 8X10	25.00					
TLUN6503	0.9 MM METRIC ALLEN KEY	25.00	X	X	X	X	X
TLAC000110	HOLDER FILM SET	164.28				X	X
TLUNAC097601	CLEANING PLATE GUIDE AY	410.10		X	X	X	X
TLUNSK000001	CALIBRATION KIT	1,066.64		X	X	X	X
TLUNSK000016	LOADER PIN HEIGHT ADJ GAUGE	690.00	X	X	X	X	X
TLUNSK000017	LOOP HEIGHT GAGE/AGFA PLATE	289.14	X	X	X	X	X

TLUNSK000024	TROLLEY-ROLLER DISTANCE GAUGE.	138.00	X	X	X	X	X	
TLUNSK000041	ROLLER SENSOR ADJUSTMENT	162.00			X	X	X	
TLUNSK000098	CR360 TOOL KIT	1,327.68						X
TLUNSK000109	AUTOLOOP SOLENOID GUAGE	96.00				X	X	
TLUNSK000152	CASSETTE EJECT POS TOOL KIT (CR360)	650.08						X
TLUNSK000153	AUTOLOOP OFF GAUGE (CR360)	346.16						X
TLUNSK000154	AUTOLOOP ON GAUGE (CR360)	346.16						X
TLUNSK250056	SCREEN EXTRACTOR	219.06	X	X	X	X	X	X
TLUNT0090252	PLATE GUIDE TOOL	142.86	X	X	X	X	X	X

Part Numbers followed by a * have a new part number, which follows in **red**

Machine	Model	M6	M7	M20	M35A	M35A-M	2000/A	MIN-R	180LP	270RA	460RA
Entrance	OLD-LARGE	458953	458953	458953	287041	287041	287041	-	458953	458953	458953
	NEW-LARGE	T990242	T990242	T990242	T287041	T287041	T287041	7C8322	T990242	T990242	T990242
Crossover	OLD-SMALL	458954	458954	458954	287045	287045	287045	287045	562389	562389	458954
	OLD-LARGE	458953	532020	532020	287049	287049	287049	287049	T5B6001	T5B6001	458953
	NEW-SMALL	T990243	T990243	T990243	T990305	T990305	T990305	T990305	T562389	T562389	T990243
	NEW-LARGE	T990242	-	-	T990308	T990308	T990308	T990308	T5B6001	T5B6001	T990310
Rack Side	OLD	458228	458228	458228	287006	-	287006	-	914703	914703	458228
	OLD	-	-	-	-	-	-	-	562765	562765	-
	NEW	T5B6409	T5B6409	T5B6409	-	-	-	-	T5B6411	T5B6411	T5B6409
	SMOOTH	-	-	-	-	T7C7860	-	3E0970	-	7C7859	-
Rack Bottom	OLD	562765	562765	562765	T1C7635	-	-	-	562765	562765	562765
	SMOOTH	-	-	-	-	T7C7860	-	-	-	-	-
	NEW	T5B6411	T5B6411	T5B6411	T1C7635	-	T1C7635	-	T5B6411	T5B6411	T5B6411
	SMOOTH	-	-	-	-	T7C7860	-	3E0970	-	-	-
	WASH RACK	-	-	-	-	-	-	-	8B7377	914703	-
	SUB-WASH	-	-	-	-	-	-	-	-	-	-
Squeegee	OLD	458954	-	-	-	-	-	-	-	-	458954
	NEW	T990243	-	-	-	-	-	-	-	-	T990243
Dryer	BOTTOM	-	-	-	-	-	-	-	914704	914704	-
	OLD	-	532149	-	-	-	-	-	-	-	-
	-	-	-	-	-	-	-	-	-	-	-
	SIDE/TOP	-	-	-	-	-	-	-	-	-	-
	NEW	-	-	T990243	T990305	T990305	T990305	T990305	-	-	-
	OLD	-	532101	458954	287045	287045	458954	287045	-	-	-
INTEGRATED	-	-	-	-	-	-	-	-	914705	-	
Exit Rack	TOP	-	-	-	-	-	-	-	911748	911748	-
	BOTTOM	-	-	-	-	-	-	-	911749	911749	-

