

Kodak DIRECTVIEW CR Software Options

Expand your KODAK DIRECTVIEW CR System capabilities with a full suite of enhanced optional features

Increase the productivity and ease of use of your KODAK DIRECTVIEW CR System with a choice of software options. These software packages enhance the standard KODAK DIRECTVIEW CR Software on your KODAK DIRECTVIEW CR System. Each package is offered individually—so you can choose and purchase only the option or options right for you and the way you work.

KODAK DIRECTVIEW EVP Software

- ▶ Take image diagnostic quality to a new level by increasing latitude while preserving the contrast of image detail.

KODAK Black Surround/Masking Software

- ▶ Provide better, easier viewing of CR images by masking the white collimated area to eliminate foreground flare.
- ▶ Automatically identify the area of interest and hide surrounding space.
- ▶ Manually select complex anatomy to provide proper collimation verification.

KODAK DICOM Modality Work List

- ▶ Connect the CR to the HIS/RIS for automatic transfer of patient and procedure information.
- ▶ Boost technologist productivity by eliminating manual data entry.

Grid Detection and Suppression Software

- ▶ Automatically detect and suppress gridline artifacts to improve image quality.

Low-Exposure Optimization Software

- ▶ Reduce the appearance of quantum noise in areas of low exposure while preserving diagnostic detail.

IHE Scheduled Workflow Software

- ▶ Automatically notify the HIS/RIS of exam status.
- ▶ Increase productivity with consistent patient and ordering information.



Administrative Analysis and Reporting Software

- ▶ Retrieve/consolidate information from all KODAK DIRECTVIEW CR Systems and export to any networked PC. Centralized information management enhances operations and performance.
- ▶ Improve effectiveness by tracking and trending cassette use, technologist statistics, exposure index, and other critical metrics.
- ▶ Track the frequency and cause of rejected images to improve overall performance.

Procedure Mapping

- ▶ Speed data entry with RIS and procedure codes.
- ▶ Save time and keystrokes by defining and cataloging frequently used procedures.

KODAK DIRECTVIEW CR Long-Length Imaging System Software

- ▶ Perform long-bone imaging procedures using upright or supine equipment.
- ▶ Fully automated stitching eliminates the need to orient, order, and align images.

Advanced Trauma

- ▶ Quickly select a custom patient record for emergencies when more than one trauma patient must be imaged.
- ▶ Set up and use predefined exam/procedure codes to speed trauma workflow.

KODAK DIRECTVIEW CR Mammography Feature*

- ▶ Present images from KODAK mammography cassettes using 50 micron pixels.

Remote Patient Data Entry Software

- ▶ Reduce radiographer workload by enabling data entry at the office PC.

Security Audit Log Software

- ▶ System login feature prevents unauthorized users from accessing data stored in the CR system, and tracks image access queries in support of privacy regulation requirements.

Consulting and training

- ▶ Extensive application consulting and training services.

**Not available in the U.S.*



Kodak

Licensed Product

Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608

Carestream Health Canada Company
6 Monogram Place
Suite 200
Toronto, Ontario M9R 0A1
CANADA

DirectView is a trademark of Carestream Health, Inc.
The Kodak trademark and trade dress are used under license from Kodak.
M1-311 Printed in U.S.A. 9/07 © Carestream Health, Inc., 2007 CAT No. 101 6468

More information

To learn more about these optional software packages for your KODAK DIRECTVIEW CR System, or to learn about other available software options, contact your Carestream Health representative, or call 1-877-865-6325, ext. 227.

→ www.carestreamhealth.com

Carestream
HEALTH

The Innovation Powering **Kodak** Health Products