

Kodak DIRECTVIEW CR Software

Powerful software maximizes the performance of your KODAK DIRECTVIEW CR System

KODAK DIRECTVIEW CR Software now includes additional workflow-enhancing features that are designed to enhance KODAK DIRECTVIEW CR System productivity, security, convenience, flexibility, and image quality.

Standard with all new KODAK DIRECTVIEW CR Systems, KODAK DIRECTVIEW CR Software is also available as an upgrade for existing systems.

Improve workflow

- ▶ **Digital markers**—Improve communications between radiographers, radiologists, and referring physicians using 48 electronic markers and “free text” comment areas.
- ▶ **Configurable text boxes**—Apply internal and external text boxes to single and multiforamt print images, enabling supplemental information to be included as part of the image.
- ▶ **Zoom-and-pan motion detection**—Increase image quality and reduce retakes. Magnify the image area to full resolution to verify proper positioning and to detect patient motion.
- ▶ **“Accept all images”**—One button enables all images in a study to be accepted without accessing each individual image.
- ▶ **True-size printing**—Print images at true size using KODAK DRYVIEW Laser Imagers.
- ▶ **Cancelled order support**—Automatically remove cancelled HIS/RIS patient records from the CR database.

Safeguard patient images and information

Guard the privacy and integrity of images and information with features that support security initiatives.

- ▶ **Cisco Security Agent (CSA)**—Improved protection from malware. CSA prevents malicious behavior originating from known and unknown security threats.
- ▶ **System access security**—Assign user names and passwords for access control. Guard against unauthorized viewing of an unattended system with a security screensaver.
- ▶ **Firewall protection**—Minimize threats from worms and other malware with a built-in firewall.





Feature-rich user interface

Support operating convenience and efficiency with user interface features designed to enhance workflow.

- ▶ Choice of 10 languages—Choose Chinese, English, French, German, Italian, Japanese, Norwegian, Portuguese, Spanish, or Swedish.
- ▶ Easy key operations transfer—Save key settings to disk, simplifying their transfer to other KODAK DIRECTVIEW CR Systems.
- ▶ Patient/exam required fields configuration—Ensure exam records are complete by selecting which data is mandatory.
- ▶ Patient data entry screen—Access laterality, procedure code, and patient information fields.
- ▶ System status screen—Check erase lamps, disk and memory usage, and database size statistics.
- ▶ Utility menu—Execute CR system shutdown and recovery, status functions, logout, and password changes.
- ▶ Extended patient and exam information—Record supplemental information including radiographic technique and patient notes.

Choice of software options

- ▶ Expand the capability of your KODAK DIRECTVIEW CR System with a variety of available software options.

Proactive service

- ▶ Preventive maintenance alerts.
- ▶ Secure remote diagnostics and service.

More information

For more information about KODAK DIRECTVIEW CR Software and other customer-focused options designed to enhance the performance of KODAK DIRECTVIEW CR Systems, contact your Carestream Health representative, or call 1-877-865-6325, ext. 227.

➔ www.carestreamhealth.com

Kodak

Licensed Product

Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608

Carestream Health Canada Company
6 Monogram Place
Suite 200
Toronto, Ontario M9R 0A1
CANADA

Carestream 
HEALTH

The Innovation Powering **Kodak** Health Products