


Kodak DIRECTVIEW Capture Link System

Centralized DR and CR workflow has never been more flexible, convenient, and productive

The KODAK DIRECTVIEW Capture Link System links up to five KODAK DIRECTVIEW DR and/or CR Systems to improve workflow and productivity in high-volume areas. Integrate DR and CR systems together to streamline workflow when a study includes DR and CR images. Share cassette identification, scanning, and image review functions on linked CR systems. Add up to ten KODAK DIRECTVIEW Remote Operations Panels (ROPs) or PCs with KODAK DIRECTVIEW Remote Patient Data Entry Software (RPDES) for even greater flexibility.



Link multiple DR and CR systems to improve workflow and productivity

Managing studies with DR and CR images has never been easier

- ▶ Link DIRECTVIEW DR and CR systems together seamlessly to streamline workflow and improve productivity
- ▶ Identify CR cassettes at the DR console and scan them at a remote or in-room CR reader
- ▶ View the entire study at the DR console—DR and CR images—at once for greater productivity and convenience
- ▶ Reprocess, approve, and redirect CR images from your DIRECTVIEW DR system operator console
- ▶ Forward DR and CR images together to the final output destinations

Harness the power of multiple systems

- ▶ Flexible configuration fits the way you work
- ▶ Support better workflow and increased throughput in mid- to high-volume, centralized radiography applications
- ▶ Modify linked systems to meet changing workflow needs—add or remove DR and CR systems, ROPs, or RPDES connections
- ▶ DIRECTVIEW DR and CR systems share a consistent, intuitive user interface to minimize training

Kodak DIRECTVIEW Capture Link System

Identify, scan, and review CR images *anywhere*

- ▶ **Save** radiographer time and increase productivity for better patient care
- ▶ **Identify** cassettes at up to five linked DIRECTVIEW CR systems—and up to ten wall-mounted DIRECTVIEW ROPs
- ▶ **Scan** cassettes quickly at any linked DIRECTVIEW CR system. Cassettes from one study can be scanned at one DIRECTVIEW CR system or at multiple systems*
- ▶ **Review** study images, check for quality, and deliver images to final destinations from any linked DIRECTVIEW CR system monitor or ROP

*Cassettes from a long-length image must be scanned on the same CR system

Improve CR data sharing and access

- ▶ Access worklists and all patient and exam information from any linked DIRECTVIEW CR system or device
- ▶ Update a single patient record, reprocess images, check for quality, and deliver to the output destinations from any linked CR system or device

KODAK DIRECTVIEW Capture Link System requirements

- ▶ KODAK DIRECTVIEW CR systems must be running KODAK DIRECTVIEW Version 5 Software on a 3 Ghz (or higher) PC
- ▶ All linked systems must use the same software options, body part and position buttons, user interface language, and system configuration
- ▶ Simple Capture Link:
 - Links two CRs or a CR and DR systems
 - With both systems running KODAK DIRECTVIEW V5 Software (or higher), one system is selected to be the server for the two units
 - Includes all the functions of KODAK DIRECTVIEW Capture Link
 - Allows addition of up to ten Remote Operations Panels

KODAK DIRECTVIEW
DR 9500 System



KODAK DIRECTVIEW
Elite CR System



Kodak
Licensed Product

Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608

Carestream Health Canada Company
6 Monogram Place
Suite 200
Toronto, Ontario M9R 0A1
CANADA

More information

Learn more about the KODAK DIRECTVIEW Capture Link System and other computed radiography systems, digital radiography systems, and services. Contact your Carestream Health representative, or call 1-877-865-6325, ext. 227.

➔ www.carestreamhealth.com

Carestream 
HEALTH