



## FireCR Veterinary CR Scanner

**A full-featured CR system at an affordable price.  
Built on pioneering technology to help increase  
your productivity while simplifying workflow.**

### High Throughput

Capable of reading up to 70 imaging plates/hour the FireCR delivers images in record time ready for viewing, enhancement, duplication, archiving or distribution.

### Compact and Lightweight

The FireCR is an extremely compact and light-weight CR reader designed to fit into any veterinary clinic. Set it on a desktop, mount it on the wall or move it directly to the animal when needed, even for off-site imaging. Experience the ease of this highly versatile system, capable of efficiently handling the everyday needs of professionals in veterinary clinics of any size.

### Optimized Imaging Parameters

The FireCR can be configured for nearly every clinical application, and is designed with full DICOM connectivity. All imaging parameters are optimized, resulting in digital images that can be enhanced, enlarged, duplicated, and sent to any location in seconds with no loss of resolution.

### Upgradability

FireCR scanners can be upgraded in the field, without disrupting the workflow. Customers simply invest in the scanner that meets current needs and budgets, and then scale up as their imaging volume increases. This unique feature provides maximum flexibility for veterinary clinics, and ensures that their technology investment is protected feature provides maximum flexibility for practices and clinics, and ensures that their technology investment is protected.

### Available in Different Models

- FireCR Veterinary-80 CR Scanner
- FireCR Veterinary-60 CR Scanner
- FireCR Veterinary-40 CR Scanner
- FireCR Veterinary-20 CR Scanner



### Highlights

- High throughput - up to 70 plates per hour
- Ultra-light - 65 lbs (30 kg) CR reader
- Ultra clean and crisp images, advanced image processing
- Field upgradeable
- Compact self-contained table top or portable system
- Wall mount option for minimum space occupancy
- One step integrated scanning and erasing
- Automatic cassette and plate handling
- Rugged, torsion stabilized optics
- Delivered in a solid, easy-to-transport wheeled flight case



# FireCR Veterinary CR Scanner with PACS Pro™ Acquire



**FireCR Veterinary CR Scanner Models**

- Veterinary-80 CR Scanner** • Up to 70 plates/hour (14x17)
- Veterinary-60 CR Scanner** • Up to 60 plates/hour (14x17)
- Veterinary-40 CR Scanner** • Up to 40 plates/hour (14x17)
- Veterinary-20 CR Scanner** • Up to 25 plates/hour (14x17)



**Vertical Wall Mount**

- Offers a footprint of less than 5 inches
- Allows installation in Veterinary clinics with limited space



**Transport Case**

- FireCR is delivered in a solid transport case on wheels
- Provides veterinary clinics with a portable solution

## FireCR Specifications

### Supported IP sizes:

- Universal Cassette 14x17: Supports IP sizes 14x17" & 14x14"
- Universal Cassette 14x10: Supports IP sizes 10x12", 8"x10", 15cm x 30cm, 20cm x 30cm, 25cm x 30cm

*\*System Includes one each: 14x17" & 10x12" Cassette/IPs*

**Resolution:** Data Capture: 16-bits per pixel, 65,000 gray tones

**Dimensions (HxWxD):** 4.6" x 18.3" x 27.7" (117 x 464 x 703 mm)

**Weight:** 65 lbs (30kg)

### System set-up options

- Desktop
- Vertical wall mount

### PACS Pro Veterinary Acquire™ Software

- Intuitive Patient Worklist GUI for Patient Entry
- Acculmaging™ Auto-Recognition Image Processing
- DICOM Store/Send to any DICOM destination
- Auto & Manual Image Routing Settings
- Measuring Tools, Annotations & Study Notes
- Single & Multi-Patient CD/DVD Burning
- Touch Screen-capable Interface
- DICOM Print, PACS & Additional Viewers available

### Power requirements

- 110/220/240V, 60/50cs, 150W
- UPS Required

### Computer workstation

- Intel Core Duo/Core 2 Processor
- 320 GB Hard Disk, 2 GB RAM
- 19" Widescreen LCD 1.3MP, DVI Interface
- Windows 7 Professional Operating System
- 16xCD/DVD+/-, NetCard (100 MBit), 2 USB 2.0 ports

### Regulatory approvals

- CE (0120)
- UL, cUL, FCC
- FDA

### Safety Standards

- EN 60601-1
- EN 60825-1
- EN 60601-2

### Operating conditions

- Temperature: 59-86 Fahrenheit (15-30 Celsius)
- Temperature Gradient: 0.5 Deg. Celsius / Min.
- Humidity: 15%-95% RH
- Magnetic Field: Max. 1260 Micro Tesla

**For Sales Please Contact Your Local X-Ray Service Provider**