



Make a cost-effective upgrade to your existing KODAK DirectView CR systems to perform both mammography and general radiography imaging.

Make the most of your investment in a KODAK DirectView CR System and capitalize on a cost-effective transition to digital mammography with the KODAK DirectView CR Mammography Feature.

This upgrade lets you add mammography imaging capability to DirectView CR 850, 950, 975, Elite, or Classic systems while retaining all the features and capabilities needed for general radiography exams—without any change in your CR procedures, workflow, or user interface.

### Advantages at a glance

- Easy transition to digital mammography
- Mammography quality images for screening and diagnostic exams
- Cassettes and screens designed for mammography
- Part of a complete mammography workflow solution

# Digital mammography with Carestream Health



## An easy transition to digital mammography

- Add the CR Mammography Feature with as a retrofit to an existing system
- An upgrade delivers full mammography capability
- Compatible with existing mammography x-ray units and workflow
- Ideal for hospitals, imaging centers and clinics which perform both mammography and general radiography imaging

## Mammography-quality images for screening and diagnostic exams

- 48.5-micron pixel pitch for high-resolution scans
- Wide exposure latitude helps reduce retakes
- Breast periphery and skinline visualization means no loss of breast parenchyma
- KODAK DirectView EVP Plus, Low Exposure Optimization, Grid Detection/Suppression, and Black Surround/Masking software provides optimal quality for mammography

## Cassettes and screens specifically designed for mammography

KODAK DirectView CR Mammography Cassettes with EHR screens feature high signal and low noise.

## Additional CR accessories for mammography imaging

- KODAK DirectView CR Mammography Total Quality Tool—Streamline your CR mammography testing without the cost, scheduling, or disruption of third-party testing.
- KODAK DirectView CR Remote Operations Panel lets you communicate with up to eight KODAK DirectView CR systems.
- KODAK DirectView Capture Link System links multiple DirectView CR systems for streamlined workflow and improved productivity in high-traffic areas.
- With over 3,000 sites worldwide using the CR Mammography Feature, you can depend on Carestream Health's experience in mammography imaging and CR technology to bring digital mammography to your facility.

## Discover a complete mammography workflow solution

- CARESTREAM Mammography Workstation—a true multi-modality solution—displays images and data from a full range of modalities, including FFDM, CR, US, MR DR, NM, and PET.
- KODAK DryView 6800 Laser Imager delivers 650 ppi resolution and printing flexibility for facilities performing both mammography and general radiography exams. Built-in compatibility with the CR Mammography Feature.
- KODAK DryView 8900 Laser Imaging System Mammography Upgrade enables 650 ppi
- KODAK CARESTREAM Information Management Solutions lets you effectively manage the exploding volumes of digital mammography images and information.



Want to find out more?

<http://mammo.carestreamhealth.com>  
1-888-777-2072

# Carestream

Integrated solutions for every step of diagnostic imaging  
**CR+DR+RIS+PACS+ARCHIVE+PRINT+**