

## MAVIG Update – IPL 2014/1 Technical Products



Protection & Suspension Systems  
Pricelist & Portfolio IPL 2014/1 TP

Dear MAVIG Customers,

Our product offering has seen quite a lot of positive changes this year. To keep our customers up to date on all the changes, we have issued a new portfolio and price list.

We are pleased to say that the new price list was not created because of a price increase, but is an expansion from the IPL2013/2 TP to include all of our new products.

Please take the time and find out about the changes and new additions to MAVIG's portfolio. Also, feel free to order a printed or digital copy of the new "Portfolio" and "Pricelist".




### Innovative Solution for Stationary Mounted Systems

#### Portegra2 TriPin Ceiling Column

The new TriPin ceiling column offers the possibility to mount three suspension arms. All of Portegra2's features and slim design have been kept, offering a solid price competitive option to other large central axis options.

See page 11



### Substructure Size Change

#### Technical Change

The range of the substructure for the Portegra2 has been increased. The minimal length has now been shorted to 200 mm.

See page 11



### New Concept for an Optimized, above Table, Scattered Radiation Protection

#### OT54001 and Protective Drapes

MAVIG's new development is a new, innovative radiation protective shield with an overlapping panel curtain and ergonomically fitted radiation protective drapes, or alternatively with our practical disposable protective drapes.

See pages 12-14

## MAVIG – Update: IPL 2014/1 Technical Products

NEW



### OT50001 with another Connection Piece Length

#### OT50001V60

Our traditional protection shield OT50001 has an additional connection bar length.

You now have the option to choose a shorter 600 mm bar instead of the standard 800 mm.

See page 13

NEW



### Protective Shield OT50001 „Mirrored“

#### OT50001I

Our well known shield design with couture cut-out on the right side for cath lab procedures now comes in a „mirrored“ version so that it can be placed on the other side of the table.

The connection bar still faces towards the user for proper positioning.

See page 13

NEW



### Contour-Fitted, Sterile Covers

#### for OT50001, OT540001, WD261, and Radiation Protective Drapes

New for our portfolio are sterile covers for our shields OT50001 and OT54001. These covers are contour-fitted and cover the complete shield.

Also new are covers for our WD261 and multiple-use patient drapes, which further optimize their use.

See pages 13, 14, and 37

NEW



### 160,000 Lux LED OR Lamp

#### LED5 SC and LED5 MC for Hybrid Central Axis Systems

For optimal light output in the OR, MAVIG offers a 160,000 Lux lamp.

Depending on the requirements, the lamp is available as a single and a multi-colour version.

See page 17

NEW



### HD and SD Camera

#### KAM01HD & KAM01SD

You now have the choice of ordering a high-resolution HD or a normal SD camera with our special precabled LED3 MC OR lamp with 130,000 Lux.

See page 17

## MAVIG – Update: IPL 2014/1 Technical Products

NEW



### New Hybrid Central Axis - Long Arm Ceiling Suspension

#### For Equipment up to 21 kg

Our new central axis system has an extreme wide reach to make sure that ceiling hung devices can be placed where they need to be.

This is especially important for rooms like a Hybrid OR.

See page 19

HZA

NEW



### Optimized Lower Body Radiation Protection

#### UT60 Revised

The UT60 is MAVIG's answer to table mounted protection for tiltable tables. The lower sway panels self-adjust to a table inclination of up to 15°.

The modified panel lengths provide optimal protection.

See pages 32 and 33

UT6001

NEW



### Additional Accessories Rail for Lower Body Protection Adapter

#### More Space

The UT60 and UT69 table mounted shields can now be ordered with an additional rail, which fits over the adapter's metal tightening knob.

This rail can now be used to hang small accessories like an additional protective panel or small joystick.

See pages 32 and 33

Additional Rail

NEW



### Body-Contoured, Height-Adjustable, Encapsulating, Mobile Radiation Shield

#### WD261

The new WD261 mobile shield developed by MAVIG provides mobile X-ray protection in any location where X-ray protection is currently absent or lacks flexibility, leaving your safety at risk and disrupting your work routine.

See pages 36 and 37

WD261

## MAVIG – Update: IPL 2014/1 Technical Products



### Halogen Lamps

#### Discontinued

LEDs have now taken over halogen lamps as the industry standard.

Due to the superiority of LEDs and minimal demand for halogen lamps, MAVIG has taken these lamps out of its program.

Halogen



### Ceiling Suspended Shield OT50010

#### No longer in our Portfolio

The protective shield OT50010 has been replaced by the new OT54001 shield, which offers better protection and can be used for both femoral and radial access interventions.

OT50010



### Protective Curtain VH192

#### Discontinued

Due to low demand, it was no longer economical to continue to manufacture this model.

The VH192 has been taken out of our active portfolio.

VH192



### Monitor Holders GD65, GD72/76 & MD70

#### Regulatory Affairs

Due to changing regulatory requirements, these products may not be listed for all countries.

If you have any questions about specific country requirements, please feel free to contact our sales department.

See pages 26, 27, and 30

Regulatory

OT50002	Transparent acrylic shield OT50002, dimmable, height adjustable handle.
OT50003	Transparent acrylic shield OT50003, dimmable, with a patient contour cut-out for positioning.
OT50011	Transparent acrylic shield OT50011, dimmable, rubber panels in hygienic PVC covers.
OT34001	Transparent acrylic shield OT34001, dimmable, cut out and curtain for medial and lateral use. Includes a preset Portegra2 suspension arm.
STRA-014	50 sterile, disposable full covers for OT50002.
OT80001	Transparent acrylic shield OT80001, dimmable, with flexible lead rubber panels in hygienic PVC covers.
OT80002	Transparent acrylic shield OT80002, dimmable, with flexible lead rubber panels in hygienic PVC covers.

### Reference of the Portegra2 Component Parts

#### Better Referencing

Certain sales articles are made up of more than one component. Each of these parts has their own serial number and/or regulatory listing.

In order to help you keep a clear reference of all of these parts, these component part numbers found on the labels have been listed.

Reference

**MAVIG GmbH**  
Headquarters

Stahlgruberring 5  
81829 Munich  
Germany

Phone +49 (0) 89 420 96 0  
Fax +49 (0) 89 420 96 200  
eMail info@mavig.com