

Carestream



Vita XE CR Systems

DIGITAL. COMPACT. AFFORDABLE.

CARESTREAM DIRECTVIEW Vita XE CR Systems are small-footprint solutions that provide high-quality digital images – at a very affordable price. With up to 69 plates per hour, Vita XE CR Systems boost efficiency and productivity, while supporting enhanced patient care.



EFFICIENCY DEFINED.



CARESTREAM Vita XE CR Systems

BIG PERFORMERS FOR SMALLER FACILITIES.

Vita XE CR Systems are designed for small-to-midsized hospitals, clinics and practices that need to watch every dollar spent – but won't compromise on the quality of their images.

THE RIGHT SIZE. THE RIGHT SOLUTION.

Why pay for capabilities you don't need and won't use? Many large, complex digital-imaging systems are overkill for your facility's applications. And, they carry a prohibitive price tag. Vita XE CR Systems provide:

- Quality imaging at an affordable cost
- The ability to offer in-house digital X-ray
- Faster turnaround time – for speedier diagnoses and enhanced care






EASY SETUP. QUALITY RESULTS.



SCHEDULE




ACQUIRE



VIEW/EDIT



DIAGNOSE/REPORT



PRINT/ARCHIVE

Set Up and Start Scanning in Just Hours

Fast, easy setup will have your system up and running the same day plus training is minimized with our intuitive interface.

Boost Your Throughput and Productivity

Pre-configured image-quality “looks” simplify selection of your default image preference plus our “smart erase” feature optimizes erase time.

Simplify Procedures and Cut Cost

Say good-bye to purchasing film and storing chemicals.

Obtain High Quality Images

Simple, easy to use image acquisition software with Enhanced Visualization – Automatic Image Processing (EVP+).

Perform Specialized Applications

LLI, Dental, Military, Mobile, Chiropractic, and Orthopedic with specially designed hardware and/or measurement tools.

HIGH-PERFORMANCE DIAGNOSTIC IMAGING.

CARESTREAM CR Screens and Cassettes

Ideal for use with the Vita XE CR Systems, our line of CR cassettes and phosphor screens offers a wide selection and high performance to meet the needs of various CR exams and clinical applications, including long-length imaging, orthopaedics, pediatrics, and more.

- 8 x 10 in.
- 10 x 12 in.
- 11 x 14 in.
- 14 x 14 in.
- 14 x 17 in.
- 14 x 33 in.
- 24 x 30 cm
- 15 x 30 cm



System Specifications

DIMENSIONS

34 cm x 75 cm x 55 cm / 29.5 in x 21.7 in x 13.4 in

WEIGHT

36 kg / 79 lbs

SCANNING RESOLUTIONS

12 Pixels/mm

6 Pixels/mm

FEATURES

Enhanced User Interface for Acquisition

Web Order Entry

DICOM and HL7 Compliant

DICOM Print

Multi-format Printing

DICOM Viewer

AVAILABLE OPTIONS

Mini-PACS Archiving

Clinical Reporting

Web Viewer License

Long-Length Imaging

DICOM Store Service Class Provider (SCP)

Advanced Measurement Tools Set

OmniLink Software

Tablet Viewer

Dental Panoramic Package

A Community of Service and Support.

Tap in to all the advantages of our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver – based on thousands of customer engagements worldwide and our 100-year heritage in medical-imaging innovation.



TI-BA Enterprises, Inc. is a master distributor for Carestream Directview Vita products. We provide sales and technical support on all Vita systems, as well as parts and bench service for dealers and end-users.

Call TI-BA at (585) 247-1212 or visit: www.ti-ba.com

TI-BA Enterprises, Inc. - 25 Hytec Circle, Rochester NY 14606

carestream.com/vita

Carestream

