

## Splash-proof DR Panel Protector

The 14×17 Splash Proof DR Panel Protector is a full enclosure that provides additional fluid and drop protection to sensitive flat panel detectors. With quick release locks, the case is both efficient and easy to use allowing the radiographer to focus on positioning and patient care instead of panel protection.

Ideal for the any environment where a panel will be handled regularly and subjected to potential fluid damage or accidental drops.



### Product Advantages

---

- Particulate and dust protection
- Drop protection
- Quick release locks
- Curved corners
- Weight Tolerance, 400lbs.
- Smooth non-porous surface
- Lightweight
- Optional cutouts for LEDs and power button
- IP5X Rating for dust per IEC 60529
- Fully encases panel in a plastic shell