



# LG Veterinary Digital X-ray Systems



17HK700G-W  
 17HK701G-W | 14HK701G-W

## Superior Image Quality for Diagnostic Confidence

Highly Sensitive Csl AED<sup>\*</sup> Detector with 16bit AD Conversion (Data Input & Output) \*Auto Exposure Detection



## Light yet Durable

**IP41: Water Resistance & Dustproof**



**Solids Level 4**  
 Protected against objects of size 1mm or larger



**Liquids Level 1**  
 Protected against dripping liquids, such as blood and bodily fluids

**Magnesium & Carbon-fiber Body**

With a super strong body that combines carbon-fiber with magnesium, The LG DXD is not only lightweight, but also exceptionally strong.

**Product Weight**

17HK700G-W    17HK701G-W    14HK701G-W  
Including Battery    Including Battery    Including Battery



**Maximum Load Weight**

17HK700G-W & 17HK701G-W  
 14HK701G-W




**Uniform load**  
 150kg (331lb)  
**Point load**  
 100kg (220lb)

## Acquisition Workstation Software (AWS) for LG DXD


**Veterinary User Interface**

LG AWS offers your practice an easy-to-learn Graphical User Interface (GUI) with animal specific exam views to make order entry and patient imaging simple and efficient. Pre-programmed views for canine, feline, equine and exotic pets.



## Specification

### Digital X-ray Detector

Model Name		17HK700G-W (Wired)	17HK701G-W (Wired / Wireless)	14HK701G-W (Wired / Wireless)
External Dimension / Weight		460mm x 460mm x 15.5mm 	460mm x 460mm x 15.5mm 	384mm x 460mm x 15.2mm 
		3.2kg, 7lb	3.4kg, 7.5lb (Incl. Battery)	2.95kg, 6.5lb (Incl. Battery)
Intended Use		General Radiography		
TFT Type		A-si		
Scintillator		CsI		
Pixel Size		140μm		
Image Area	X-ray Sensitive Pixel	3072 x 3072 pixels	3072 x 3072 pixels	2500 x 3052 pixels
	X-ray Sensitive Area	430.08 x 430.08 mm	430.08 x 430.08 mm	350 x 427.28 mm
MTF	MTF @ 0.5lp/mm	(Typ.) 89%		
DQE	DQE @ 0.1lp/mm	(Typ.) 72%		
A/D Conversion		16bit		
Data Output		16bit		
Communication	Type	Wired	Wired / Wireless	Wired / Wireless
	Wireless	-	IEEE802.11a/b/g/n/ac, 2.4GHz / 5GHz	IEEE802.11a/b/g/n/ac, 2.4GHz / 5GHz
	AED (Auto Exposure Detection)	Support		
Cycle Time		8sec.	8sec. (Wired) / 11sec. (Wireless)	8sec. (Wired) / 11sec. (Wireless)
Full Image		3sec.	3sec. (Wired) / 6sec. (Wireless)	3sec. (Wired) / 6sec. (Wireless)
Durability	Maximum Load Weight	Full (Uniform Load) : 150kg, 331lb Local (Point Load) : 100kg, 220lb		
	Resistance Rating	IP41		
	Sensor Protection Plate	Carbon Fiber Plate		
Battery	Operating Time	-	(Typ.) 240 shots / 6.0hrs. (Condition : Interval Time 90sec.)	(Typ.) 240 shots / 6.0hrs. (Condition : Interval Time 90sec.)
	Duration in Standby Status	-	(Typ.) 13hrs.	(Typ.) 13hrs.
	Charging Time	-	3hrs.	3hrs.

### Acquisition Workstation Software

Installation Requirement	OS	Microsoft Windows 7, 10(64 bit)
	DB	MariaDB
	Processor Options	Intel CPU Desktop 6th gen i7 or above
	RAM	8GB
	HDD	2TB

**LG Electronics Inc.**

<http://www.lge.com>

Copyright © 2018 LG Electronics. All Rights Reserved.

Distributed by

**TI-BA ENTERPRISES, INC.**

TI-BA Enterprises, Inc.  
25 Hytec Circle, Rochester, NY 14606  
Tel (585) 247-1212  
[www.ti-ba.com](http://www.ti-ba.com)  
[sales@ti-ba.com](mailto:sales@ti-ba.com)