

VXR Veterinary System

Digital Radiography System for Veterinary Practices




Transform Your Veterinary Practice with Cutting-Edge Technology

The VXR Veterinary System from DRGEM is a revolutionary digital radiography system designed specifically for veterinary practices. With the VXR system, you can say goodbye to the darkroom, processor chemicals, and film. Embrace instant, high-quality digital images of companion animals that save time, reduce labor costs, and minimize unnecessary patient dose through its fast kV rising time feature.

Key Features:

- **High-Quality Imaging:** Instant digital images with superior quality for accurate diagnosis.
- **Digital Enhancement:** Easily enhance and annotate images, providing on-screen measurements and enabling seamless image transmission for referrals and consultations.
- **Comprehensive Package:** Includes Digital Imaging Receptor, X-Ray Table, computer workstation, and a full software package for an all-in-one solution.
- **Generator Options:** Choose from a variety of generator options tailored to your practice's needs.
 - Standard: 20 / 25 / 32 / 40 kW
 - Capacitor: 20 / 25 / 32 / 40 / 50 kW



Advanced Features for Optimal Performance:

- **Movable 4-Way Tabletop:** Offers flexibility and ease of positioning.
- **X-Ray Tube Toshiba E7242X:** With a focal spot size of 0.6/1.5mm and max 200KHU anode heat storage for precise imaging.
- **GXR High-Frequency X-Ray Generator:** Ensures excellent accuracy, reproducibility, and long-term stability.
- **Comprehensive Imaging Software:** Optional Mini-PACS server & viewers for advanced image management.
- **Efficient Animal Handling:** Buckle fastener for secure animal positioning and automatic APR selection for ease of use.

Why Choose the VXR Veterinary System?

The VXR Veterinary System is more than just an X-ray machine; it's a comprehensive solution that enhances your diagnostic capabilities, improves patient care, and streamlines your workflow. With its advanced features, customizable options, and ease of use, the VXR system is the perfect addition to any veterinary practice aiming to leverage the power of digital radiography.

 ENTERPRISES, INC.

Distributed by: TI-BA Enterprises, Inc. – 25 Hytec Circle, Rochester, NY 14606
www.ti-ba.com - Tel (585) 247-1212 - info@ti-ba.com